

## SLC library pooled electroporation, plating, determination of transformation efficiency and maxi prep

### Protocol adapted from the GeCKO library amplification protocol (Zhang's lab)

1. Dilute the library to 50 ng/uL in water or TE (if not already diluted).
2. Electroporate the library
  - a. Add 2 uL of 50 ng/uL library to 25 uL of electrocompetent cells with an efficiency of  $\geq 10^9$  cfu/ug. We generally use Lucigen Endura cells (cat # 60242)
  - b. Electroporate using the manufacturer's suggested parameters/protocol.
  - c. Recover in 975 uL recovery media (i.e. media provided with cells) and transfer to a loosely capped tube with an additional 1 mL of recovery media.
  - d. Repeat for a total of 4 electroporations and rotate at 250 rpm for 1 hour at 37 C
3. Plate a dilution to calculate transformation efficiency
  - a. Pool all 8 mL of electroporated cells. Mix well.
  - b. Remove 10 uL and add to 1 mL of recovery media, mix well, and plate 20uL onto a pre-warmed 10cm petri dish (ampicillin). This is a 40,000-fold dilution of the full transformation and will enable you to estimate transformation efficiency to ensure that full library representation is preserved.
4. Plate the transformations

*Follow Step a) if your lab has 24.5 cm<sup>2</sup> bioassay plates for large-scale bacteria culture; otherwise follow Step b), which substitutes 20 standard (10 cm round) petri dishes.*

  - a. Plate 4mL of transformation on each of 2 pre-warmed 24.5 cm<sup>2</sup> bioassay plates (ampicillin) using a spreader. Spread the liquid culture until it is largely absorbed into the agar and won't drip when turned upside down.
  - b. Alternatively, spread 400 uL of transformation mix per petri onto 20 pre-warmed petri dishes.
5. Grow all plates inverted for **14 hours at 30-32 C**. Growth at this lower temperature reduces recombination between the lentiviral long-terminal repeats. (Growth at 37 C is also acceptable if 32 C is not possible.)
6. Calculate transformation efficiency
  - a. Count the number of colonies on the dilution plate.
  - b. Multiple this number of colonies by 40,000 for the total number of colonies on all plates.
  - c. Proceed if the total number of colonies is at least  $3 \times 10^6$ . This efficiency is equivalent to ~1000X colonies per construct.
7. Harvest colonies
  - a. Pipette 10 mL of LB onto each 24.5 cm<sup>2</sup> bioassay plate (or, 500 uL per 10 cm petri dish)
  - b. Scrape the colonies off with a cell spreader/scrapper.
  - c. Pipette off the liquid plus scraped colonies into a tube and repeat the procedure a second time on the same plate with additional 5-10 mL. *Note: Weigh this tube prior to adding any liquid to it.*
8. Weigh the bacterial pellet to determine the proper number of maxiprep columns to use
  - a. Spin down all liquid to pellet the bacteria and then discard the supernatant.
  - b. Weigh the bacterial pellet and subtract the weight of the tube
9. Maxi-prep for downstream virus production and future amplification
  - a. Using a maxi scale plasmid prep, each column can handle approximately 0.45 g of bacterial pellet.
  - b. Perform a sufficient number of maxi preps so as to not overload a column.