

TF ORF barcode library PCR for NGS

Introduction

PCR protocol for amplification of TF ORF barcodes from plasmid or gDNA library for NGS. For additional details on using the TF ORF library for screening, please refer to Joung et al, "Genome-scale CRISPR-Cas9 knockout and transcriptional activation screening," Nature Protocols 2017.

Materials

› Forward Primers

› ORF-NGS-F1:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTTAAAGTAGAGTTGGCTTTATATATCTTGTG
GAAAGGACGA

› ORF-NGS-F2:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTATCATGCTTATTGGCTTTATATATCTTGTG
GAAAGGACGA

› ORF-NGS-F3:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTGATGCACATCTTTGGCTTTATATATCTTGT
GGAAAGGACGA

› ORF-NGS-F4:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTCGATTGCTCGACTTGGCTTTATATATCTT
GTGGAAAGGACGA

› ORF-NGS-F5:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTTCGATAGCAATTCTTGGCTTTATATATCTT
GTGGAAAGGACGA

› ORF-NGS-F6:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTATCGATAGTTGCTTTTGGCTTTATATATCT
TGTGGAAAGGACGA

› ORF-NGS-F7:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTGATCGATCCAGTTAGTTGGCTTTATATAT
CTTGTGGAAAGGACGA

› ORF-NGS-F8:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTCGATCGATTTGAGCCTTTGGCTTTATATA
TCTTGTGGAAAGGACGA

› ORF-NGS-F9:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTACGATCGATACACGATCTTGGCTTTATAT
ATCTTGTGGAAAGGACGA

› ORF-NGS-F10:

AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCTTACGATCGATGGTCCAGATTGGCTTTATA
TATCTTGTGGAAAGGACGA

› Reverse Primers (only need 1 primer per condition, barcodes are bolded)

› ORF-NGS-R1:

CAAGCAGAAGACGGCATAACGAGAT**TCGCCTTGGT**GACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACAT
AGCGT

- › ORF-NGS-R2:
CAAGCAGAAGACGGCATAACGAGAT**ATAGCGT**CGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R3:
CAAGCAGAAGACGGCATAACGAGAT**GAAGAAGT**GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R4:
CAAGCAGAAGACGGCATAACGAGAT**ATTCTAGGGT**GACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R5:
CAAGCAGAAGACGGCATAACGAGAT**CGTTACCA**GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R6:
CAAGCAGAAGACGGCATAACGAGAT**GTCTGATGGT**GACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R7:
CAAGCAGAAGACGGCATAACGAGAT**TTACGCAC**GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R8:
CAAGCAGAAGACGGCATAACGAGAT**TTGAATAGG**TGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R9:
CAAGCAGAAGACGGCATAACGAGAT**TCCTTG**GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R10:
CAAGCAGAAGACGGCATAACGAGAT**ACAGGTAT**GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R11:
CAAGCAGAAGACGGCATAACGAGAT**AGGTAAGG**GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
 - › ORF-NGS-R12:
CAAGCAGAAGACGGCATAACGAGAT**AACAATGG**GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTTAAAGCAGCGTATCCACATAGCGT
- › NEBNext High Fidelity PCR Master Mix, 2× (New England BioLabs, cat. no. M0541)

Procedure

Library PCR for NGS

1. We have provided NGS primers that amplify the TF ORF barcode region with Illumina adaptor sequences. Pool the 10 ORF-NGS-F primers at equimolar ratio to produce ORF-NGS-F (pooled). These primers stagger the resulting PCR amplicons and diversify the library.

2. To prepare the TF ORF barcode library for NGS, set up a reaction using the ORF-NGS-F (pooled) and a unique ORF-NGS-R for each condition according to the following table. For amplifying from plasmid library, use 20ng of plasmid per reaction as template. For amplifying from gDNA, use 1ug-3ug of gDNA per reaction as template. Use a different ORF-NGS-R primer for each condition, and scale up the number of reactions according to the scale of the screen.

NEBNext PCR reaction components				
	A	B	C	D
1	Component	Starting concentration	Amount per reaction (uL)	Final concentration
2	NEBNext PCR Master Mix	2X	25	1X
3	ORF-NGS-F (pooled)	10 uM	1.25	0.25 uM
4	ORF-NGS-R (unique)	10 uM	1.25	0.25 uM
5	Template	20 ng/uL	1	0.4 ng/uL
6	H2O		21.5	
7			50	

3. Perform a PCR using the following cycling conditions:

NEBNext PCR Cycle Conditions				
	A	B	C	D
1	Step	Temperature	Duration	Cycles
2	Initial denaturation	98 °C	3 min	1
3	Denaturation	98 °C	10 sec	22
4	Annealing	60 °C	10 sec	
5	Extension	72 °C	25 s	
6	Final extension	72 °C	2 min	1

4. After the reaction is complete, proceed according to the "Next-generation sequencing of the amplified sgRNA library to determine sgRNA distribution" section in Joung, Nature Protocols 2017.