

TF ORF barcode dial out PCR (10X)

Introduction

Get started by giving your protocol a name and editing this introduction.

Materials

› Forward primers

› ORF-perturb-NGS-F: AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT

› Reverse primers

› ORF-perturb-NGS-R1:

CAAGCAGAAGACGGCATACGAGATTCGCCTTGGTGAAGTTCAGACGTGTGCTCTTCCGATCTTAAGTAGAGTTGGCTTTATA
TATCTTGTGGAAAGGACGA

› ORF-perturb-NGS-R2:

CAAGCAGAAGACGGCATACGAGATATAGCGTCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTATCATGCTTATTGGCTTTATA
TATCTTGTGGAAAGGACGA

› ORF-perturb-NGS-R3:

CAAGCAGAAGACGGCATACGAGATGAAGAAGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTGATGCACATCTTTGGCTTT
ATATATCTTGTGGAAAGGACGA

› ORF-perturb-NGS-R4:

CAAGCAGAAGACGGCATACGAGATATTCTAGGGTGAAGTTCAGACGTGTGCTCTTCCGATCTCGATTGCTCGACTTGGCTT
TATATATCTTGTGGAAAGGACGA

› ORF-perturb-NGS-R5:

CAAGCAGAAGACGGCATACGAGATCGTTACCAGTGAAGTTCAGACGTGTGCTCTTCCGATCTTCGATAGCAATTCTTGGCTT
TATATATCTTGTGGAAAGGACGA

› ORF-perturb-NGS-R6:

CAAGCAGAAGACGGCATACGAGATGTCTGATGGTGAAGTTCAGACGTGTGCTCTTCCGATCTATCGATAGTTGCTTTGGCTT
TTATATATCTTGTGGAAAGGACGA

› ORF-perturb-NGS-R7:

CAAGCAGAAGACGGCATACGAGATTTACGCACGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTGATCGATCCAGTTAGTTGG
CTTTATATATCTTGTGGAAAGGACGA

› ORF-perturb-NGS-R8:

CAAGCAGAAGACGGCATACGAGATTTGAATAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTCGATCGATTTGAGCCTTTG
GCTTTATATATCTTGTGGAAAGGACGA

› NEBNext High Fidelity PCR Master Mix, 2× (New England BioLabs, cat. no. M0541)

Procedure

Library PCR for NGS

1. We have provided NGS primers that amplify the TF ORF barcode region with Illumina adaptor sequences. Pool the 10 ORF-perturb-NGS-R primers at equimolar ratio to produce ORF-perturb-NGS-R (pooled). These primers stagger the resulting PCR amplicons and diversify the library.
2. To prepare the TF ORF barcode library for NGS, set up a reaction using the ORF-perturb-NGS-F and ORF-perturb-NGS-R (pooled) for each 10X channel according to the following table.

NEBNext PCR reaction components				
	A	B	C	D
1	Component	Starting concentration	Amount per reaction (uL)	Final concentration
2	NEBNext PCR Master Mix	2X	25	1X
3	ORF-perturb-NGS-F	10 uM	1.25	0.25 uM
4	ORF-perturb-NGS-R (pooled)	10 uM	1.25	0.25 uM
5	WTA cDNA template	10-20 ng/uL	2	0.4-0.8 ng/uL
6	H2O		20.5	
7			50	

3. Perform a PCR using the following cycling conditions:

NEBNext PCR Cycle Conditions				
	A	B	C	D
1	Step	Temperature	Duration	Cycles
2	Initial denaturation	98 °C	3 min	1
3	Denaturation	98 °C	10 sec	22
4	Annealing	60 °C	10 sec	
5	Extension	72 °C	40 s	
6	Final extension	72 °C	2 min	1

4. After the reaction is complete, proceed according to the "Next-generation sequencing of the amplified sgRNA library to determine sgRNA distribution" section in Joung, Nature Protocols 2017.