

PCR amplification

- For a 1000X coverage, 40 pg of pooled library are required.
- Primers are available in a separate Excel file (CRISPR_Libraries_NGS_primers_3-17.xlsx)

PCR reaction:

gDNA	x μ L to 40 pg
Q5	0.75 μ L
5x GC Enhancer Buffer	10 μ L
5x Reaction Buffer	10 μ L
Nucleotides 10 mM	1 μ L
10 μ M Primer F	2 μ L
10 μ M Primer R	2 μ L
H2O	x μ L (to 50 μ L tot)

PCR program:

	Temp	Time	Cycles
INITIAL DENATURATION	98°C	1min	1
DENATURATION	98°C	20s	24
ANNEALING	60°C	30s	
EXTENSION	72°C	30s	
FINAL EXTENSION	72°C	5min	1
HOLD	4°C	forever	

- If multiplexing, load 10 μ L of the (PCR purified) samples in a 1.5-2% gel for the quantification and the normalization of the amount of DNA present in the band by ImageJ (band at 280bp)
- Pool together all the samples, normalizing the amount of DNA added so that each sample has a similar amount of the ~280bp band
- Load on a 2% agarose. Cut the bands at 280bp of the gel and purify with the Qiagen kit:
 - we do not heat up the sample to melt the gel but use 100 μ L more of the required Buffer and leave it shaking, vortexing and pipetting from time to time until you get a liquid but still dense consistency
 - load the volume in 2 columns and elute with 30 μ L of H₂O/column.
- Submit 10 μ L of the purified sample (at least 0.5 ng/ μ L) for sequencing. A minimum of 64bp + 2x 8bp are required for sequencing.
- Store the remaining sample at -20C for future reference