

iChIP / enChIP to purify genomic DNA

Originally developed by Toshitsugu Fujita on September 18, 2012

Modified by Hodaka Fujii on December 29, 2013

1. Crosslinking of cells

- (1) Culture target cells. Use 2×10^7 cells (e.g. Ba/F3, DT40) for chromatin preparation.
- (2) Add 37% formaldehyde to 1% final concentration into the culture medium with cells. Incubate at 37°C for 5-10 min (usually 5 min).

Cell volume	30 ml
37% formaldehyde	810 μ l

- (3) Stop crosslinking by adding 1.25 M Glycine solution to 127 mM final concentration. Incubate at room temperature for 10 min.

Cell volume	30 ml
1.25 M Glycine	3.05 ml
1.25 M Glycine	Glycine MW: 75.07
Glycine	18.8 g / 200 ml

- (4) Collect cells by centrifugation (1,300 rpm, 4°C for 5 min).
- (5) PBS wash twice. Collect the pellet (cells). The cells can be stored at -80°C.

2. Preparation of chromatin

- (1) Suspend the fixed cells in 10 ml of CLB. Incubate on ice for 10 min.

Cell Lysis Buffer (CLB) 10 mM Tris-HCl, pH 8.0, 1 mM EDTA, 0.5% IGEPAL CA-630, 1 x protease inhibitors

	<i>40 ml</i>
1 M Tris-HCl (pH 8.0)	400 μ l
0.5 M EDTA	80 μ l
IGEPAL CA-630	200 μ l
Complete-Mini	4 tablets
DDW	39.32 ml

- (2) 2,000 rpm, 4°C for 8 min. Discard carefully the supernatant.
- (3) Suspend the pellet in 10 ml of NLB. Incubate on ice for 10 min. Vortex every 2-3 min.

Nuclear Lysis Buffer (NLB) 10 mM Tris-HCl, pH 8.0, 1 mM EDTA, 0.5 M NaCl, 1% Triton X-100, 0.5% sodium deoxycholate, 0.5% lauroylsarcosine sodium salt, 1 x protease inhibitors

	<i>40 ml</i>
1 M Tris-HCl (pH 8.0)	400 µl
0.5 M EDTA	80 µl
5 M NaCl	4 ml
Triton X-100	400 µl
10% sodium deoxycholate	2 ml
30% lauroylsarcosine sodium salt	666 µl
Complete-Mini	4 tablets
DDW	32.46 ml

10% sodium deoxycholate
sodium deoxycholate 1 g / 10 ml

(4) 2,000 rpm, 4°C for 8 min. Discard carefully the supernatant.

(5) Suspend the pellet in 10 ml of PBS. 2,000 rpm, 4°C for 10 min. Collect the pellet as the chromatin fraction. The chromatin fraction can be stored at -80°C after immediate freezing in liquid nitrogen.

3. Sonication of chromatin

(1) Suspend the collected chromatin fraction in 800 µl of MLB3. Transfer the suspension into a 1.5 ml microtube.

Modified Lysis Buffer 3 (MLB3) 10 mM Tris-HCl, pH 8.0, 1 mM EDTA, 0.5 mM EGTA, 150 mM NaCl, 0.1% sodium deoxycholate, 0.1% SDS, 1 x protease inhibitors

	<i>10 ml</i>
1 M Tris-HCl (pH 8.0)	100 µl
0.5 M EDTA	20 µl
0.1 M EGTA	50 µl
5 M NaCl	300 µl
10% sodium deoxycholate	100 µl
10% SDS	100 µl
Complete-Mini	1 tablet
DDW	9.33 ml

(2) Sonication of the chromatin by using Ultrasonic disruptor UD-201 (TOMY SEIKO).

Condition is as follows:

Output: 3

Duty: 100% (continuous)

Time: Free

10 - 18 cycles of sonication for 10 sec and cooling on ice for 20 sec

2 min on ice after 5 - 6 cycles

Keep the position of the tip of the sonication bar approximately 0.5 cm away from the tube bottom.

(3) 13,000 rpm, 4°C for 10 min. Collect the supernatant (800 µl). The supernatant can be stored at -80°C after immediate freezing in liquid nitrogen.

4. Reverse crosslinking (Evaluation of fragmentation of chromatin)

(1) Suspend 10 µl of the fragmented chromatin in 85 µl of distilled water.

(2) Add 4 µl of 5M NaCl. Incubate at 65°C overnight.

(3) Add 1 µl of 10 mg/ml RNase A. Incubate at 37°C for 45 min.

(4) Add 2 µl of 0.5M EDTA (pH 8.0), 4 µl of 1M Tris-HCl (pH 6.8), and 1 µl of Proteinase K (Roche). Incubate at 45°C for 1.5 h.

(5) Pick up 10 µl for electrophoresis in 1% agarose gel without staining dye. 100 V for 30 min.

(6) Gel staining with staining dye for 0.5-1 h.

5. Preparation of Dynabeads conjugated with antibody

(1) Transfer 30 µl Dynabeads-Protein G (Invitrogen) in a new 1.5 ml tube.

(2) Put the tube on a magnet stand and wait for 2 min. Discard the supernatant by pipetting.

(3) Add 1 ml PBS with 0.01% Tween-20. Put the tube on a magnet stand and wait for 2 min. Discard the supernatant by pipetting.

PBS	10 ml
10% Tween-20	10 µl

(4) Repeat the step (3).

(5) Add 300 µl PBS with 0.01% Tween-20 and 0.1% BSA.

PBS	10 ml
10% Tween-20	10 µl
7.5% BSA	133 µl

- (6) Add 3 µg antibody (e.g. anti-FLAG antibody Sigma F1804, control IgG). Rotate at 4°C overnight.
- (7) Spin down briefly. Put the tube on a magnet stand and wait for 2 min. Discard the supernatant by pipetting.
- (8) Add 300 µl PBS with 0.01% Tween-20. Invert several times and spin down briefly. Put the tube on a magnet stand and wait for 2 min. Discard the supernatant by pipetting.
- (9) Repeat the step (8), twice. The Dynabeads are ready for the next step.

6. Chromatin immunoprecipitation

- (1) Transfer 160 µl of the fragmented chromatin, which corresponds to chromatin extracted from 4×10^6 cells, into a new 1.5 ml tube.
- (2) Add 340 µl of MLB3 1.47% Triton X-100 (final 1%).

MLB3	1 ml
Triton X-100	14.7 µl

- (3) Transfer all (500 µl) of the chromatin solution into the tube, in which the Dynabeads conjugated with control IgG were prepared at the step 5-(9). Rotate 4 °C for 1h.
- (4) Put the tube on a magnet stand and wait for 2 min.
- (5) Transfer the supernatant into the tube, in which the Dynabeads conjugated with specific antibody (e.g. FLAG antibody) were prepared at the step 5-(9). Rotate 4 °C overnight.
- (6) Put the tube on a magnet stand and wait for 2 min. Discard the supernatant by pipetting.
- (7) Wash 1: Add 1 ml of LSB. Rotate 4 °C for 10 min. Put the tube on a magnet stand and wait for 2 min. Discard the supernatant by pipetting. Repeat wash with LSB.

Low Salt Buffer (LSB) 20 mM Tris-HCl, pH 8.0, 2 mM EDTA, 150 mM NaCl, 1% Triton X-100, 0.1% SDS

	<i>20 ml</i>
1 M Tris-HCl (pH 8.0)	400 µl
0.5 M EDTA	80 µl
5 M NaCl	600 µl
Triton X-100	200 µl
10% SDS	200 µl
DDW	18.52 ml

- (8) Wash 2: Repeat the step (7) with HSB x 2.

High Salt Buffer (HSB) 20 mM Tris-HCl, pH 8.0, 2 mM EDTA, 500 mM NaCl, 1% Triton X-100, 0.1% SDS

	<i>20 ml</i>
1 M Tris-HCl (pH 8.0)	400 μ l
0.5 M EDTA	80 μ l
5 M NaCl	2 ml
Triton X-100	200 μ l
10% SDS	200 μ l
DDW	17.12 ml

(9) Wash 3: Repeat the step (7) with LiCl buffer x 2.

LiCl Buffer 10 mM Tris-HCl, pH 8.0, 1 mM EDTA, 0.25 M LiCl, 0.5% IGEPAL CA-630, 0.5% sodium deoxycholate

	<i>20 ml</i>
1 M Tris-HCl (pH 8.0)	200 μ l
0.5 M EDTA	40 μ l
8 M LiCl	625 μ l
IGEPAL CA-630	100 μ l
10% sodium deoxycholate	1 ml
DDW	18.035 ml

(10) Wash 7 and 8: Repeat the step (7) with **TBS with 0.1% IGEPAL-CA-630** x 2.

TBS (50 mM Tris-HCl (pH 7.5), 150 mM NaCl) with 0.1% IGEPAL CA-630

(11) Elution: Add 150 μ l of 500 μ g/ml 3xFLAG peptide (Sigma, F4799) in TBS with 0.1% IGEPAL CA-630. Incubate at 37 °C for 20 min. Put the tube on a magnet stand and wait for 3 min.

3xFLAG peptide (5 mg/ml)	50 μ l	66 μl
TBS w/ 0.1% IGEPAL CA-630	450 μ l	594 μl

(12) Repeat the elution step. Total 300 μ l.

(13) Add

5 M NaCl	12 μ l	final	200 mM
0.5 M EDTA	0.6 μ l		1 mM

(14) 65°C, overnight.

(15) Add **3 μ l** of 10 mg/ml RNase A. Incubate 37°C for 1 hr.

(16) Add

10% SDS **16 μ l**

Proteinase K **10 μ l**

(17) Incubate at 45°C for 2 hr.

(18) Purify DNA using ChIP DNA Clean & Concentrator kit (Zymo Research, D5205). Add 1.5 ml of ChIP DNA Binding Buffer.

(19) Transfer mixture into a Zymo-Spin Column in a Collection Tube.

(20) 15,000 rpm, 30 sec.

(21) Discard flow-through.

(22) Wash with 200 μ l of Wash Buffer x 2.

(23) Elute DNA with 50 μ l of Elution Buffer.