

## FREQ-Seq condensed protocol

### *Locus-specific primer design*

**Forward primer:** 5'-GTAAACGACGGCCAGT + ~20 bp primer sequence homologous to region adjacent to mutation

**Reverse primer:** 5'- CAAGCAGAAGACGGCATAACGAGCTCTTCCGATCT + ~20 bp primer sequence homologous to region ~250 bp downstream of the mutation

### *Stage 1: Locus amplification*

1. Prepare crude genomic DNA extraction from sample population (boil prep or other lysis preparation)
2. Desired loci are amplified in a 50  $\mu$ L reaction using the locus-specific FREQ-Seq primers with conditions:
  - a. 5-10  $\mu$ L of crude DNA template from Step 1.
  - b. 0.5  $\mu$ M of each locus specific primer, 200  $\mu$ M dNTP, and appropriate buffering conditions
  - c. PCR thermal cycling: 10-15x cycles using conditions optimized for the locus specific FREQ-Seq primer set.
3. Primers from this first reaction are removed by:
  - a. Treatment with 5U of exonuclease I (*Exo*I) and 1U of *Dpn*I for 30 minutes at 37 °C and 20 minutes at 80 °C to heat inactivate these enzymes.
  - b. Purification using Ampure XP™ magnetic beads, eluting in 27  $\mu$ L of ddH<sub>2</sub>O

### *Stage 2: Bar-coding and library purification*

1. For amplified loci, samples can be bar-coded as follows:
  - a. For samples treated as in Step 1.3a, 25  $\mu$ L of the *Exo*I/*Dpn*I treated reaction can be mixed with 25  $\mu$ L of fresh PCR reaction mixture at 1X concentration and enrichment primers at 2X concentrations (e.g. 0.1  $\mu$ M Primer A and Primer B).
  - b. For samples treated as Step 1.3b, 25  $\mu$ L of the purified reaction can be mixed with 25  $\mu$ L of fresh PCR master mix at 2X concentration and enrichment primers at 2X concentrations (e.g. 0.1  $\mu$ M Primer A and Primer B).
2. To each reaction, add 1/10<sup>th</sup> the molar equivalent of a single FREQ-Seq bridging primer to a single reaction (ex. for ~50  $\mu$ L reaction add 10-20 ng of dsDNA FREQ-Seq primer).
3. PCR amplify the mixture under the following conditions:
  - a. Cycling conditions:
  - b. 98 °C – 60 s
  - c. 98 °C – 20 s
  - d. 53 °C – 20 s
  - e. 72 °C – 30 s
  - f. Cycle c-e 10-15X depending on the amplicon
  - g. 4 °C – <3 hrs
4. The resulting reactions can be treated with QIAGEN Buffer PB (5M Gu-HCl (aq), 30% Isopropyl alcohol) to stop polymerase activity.
5. Pool samples according to the loci of interest, followed by purification using QIAGEN columns or generic silica gel equivalents. This is a complete FREQ-Seq library.
6. (Optional) To ensure FREQ-Seq product size and purity, products can be run on a 2% agarose/1X TAE gel followed by standard gel extraction.
7. Determine the concentration of each allele pool.
8. Dilute and mix all allele pools together to give a final concentration of 10-20 nM depending on sequencing facility requirements.