

TCRAFT vector amplification

Vector pools (1A, 1B, 2A, 2B)

Contents of vector pools:

- **1A**, containing **TRBV-TRAC** vectors for TRBV2 to 7-8 (24 vectors) (chloramphenicol)
- **1B**, containing **TRBV-TRAC** vectors for TRBV7-9 to 30 (24 vectors) (chloramphenicol)
- **2A**, containing **TRBC-TRAV** vectors for TRAV1-1 to 16 (23 vectors) (kanamycin)
- **2B**, containing **TRBC-TRAV** vectors for TRAV17 to 41 (22 vectors) (kanamycin)

The following is a protocol to amplify pools with electrocompetent cells.

Materials

1. 20 ng of each library (1A, 1B, 2A, 2B)
2. ElectroMAX DH10B *E. coli* (Thermo Fisher, 18290015) (includes S.O.C. medium)
3. Electroporation cuvettes (1 mm gap) (VWR, 89047-206)
4. LB media
5. LB plates:
 - chloramphenicol (Millipore Sigma, C0378)
 - kanamycin (RPI, 25389-94-0)
6. BTX[®] ECM[®] 630 electroporator
7. 30C incubator for plates
8. 30C shaking incubator

Protocol

1. Thaw one 100 μ L aliquot of ElectroMax DH10B *E. coli* on ice for 15 minutes.
2. Pre-cool four Eppendorf tubes and four 1-mm cuvettes on ice.
3. Gently transfer 25 μ L of DH10B bacteria to each Eppendorf tube. Add 20 ng of each library to the Eppendorf tubes (one each). Mix by tapping gently.
4. Incubate on ice for 10 minutes.
5. Transfer cell/DNA mixture into a chilled cuvette and put back on ice. Gently tap the cuvette to ensure complete contact of cells across the bottom of the chamber. Avoid bubbles.
6. Electroporate using the BTX[®] ECM[®] 630 electroporator with the settings:
 - 2000 V
 - 200 ohms
 - 25 μ F
7. Immediately add 1000 μ L of S.O.C. media right after shocking and transfer all to 15 mL Falcon tubes.
8. Incubate cells for 1h at 37C in a shaking incubator (225 rpm).
9. Transfer 1 μ L to 100 μ L of S.O.C. medium and plate (LB-Carb for 1A/1B or LB-Kan for 2A/2B). Incubate overnight at 30C. This represents a 1:1000 dilution of the electroporation product.
10. Transfer the remainder to 200 mL of LB-Carb or LB-Kan and incubate for **6 to 8 hours at 30C**. Do not overgrow (avoids plasmid multimerization and recombination). Overnight cultures are typically fine but may alter or worsen the distribution of variable regions.
11. Next day: examine the plates for colonies. >25 colonies indicates >1000x coverage. An efficient transformation should generate far more than this. Continue with midiprep/maxiprep if electroporation successful (i.e. >25 colonies per plate).

pHIV-LacZ-Sapl (destination vector)


The pHIV-LacZ-Sapl destination vector can be propagated with any standard heat shock transformation protocol. This plasmid is carbenicillin/ampicillin-resistant.

Distribution analysis of TCRAFT vector pools


1. No prep: Submit vector pools to Plasmidsaurus or another whole-plasmid sequencing service.
2. Results:
 - Download raw fastq files.
 - Check gel/histogram.
3. Process raw fastq files to examine the distribution of plasmids in each pool using accompanying code available on Addgene or here: https://github.com/birnbaumlab/Gaglione-et-al-2025/tree/main/Vector_Pool_Scripts
 - a. Code outputs table and depiction (below) of read counts for each library member.
4. Ensure all plasmids are present.

Sample data:


1A (TRBV2 to 7-8):




1B (TRBV7-9 to 30):



2A (TRAV1-1 to 16):



2B (TRAV17 to 41):



Alternatively, a conventional NGS-based approach can be used by PCR-amplifying each pool and sequencing. The following primers can be used to add partial Illumina adapters.

TRBV pools 1A and 1B:

TRBV_f ACACTCTTTCCCTACACGACGCTCTTCCGATCTNNNNNNATCTAGAGCCACCGGCATGAG
TRBV_r GACTGGAGTTCAGACGTGTGCTCTTCCGATCTNNNNNNcagactgtcactggatttagagtctctca

TRAV pools 2A and 2B:

TRAV_f ACACTCTTTCCCTACACGACGCTCTTCCGATCTNNNNNNGACGTTGAGGAAAACCCAGGAC
TRAV_r GACTGGAGTTCAGACGTGTGCTCTTCCGATCTNNNNNNCACTGAGCCTCCACCTAGCC