



BD Biosciences
Clontech
Discovery Labware
Immunocytometry Systems
Pharmingen

Search the Site

GO

CONTACT US

SITEMAP

Home About eCatalog Services Support How to Order

Clontech > Techno > Vectors > VectorsE > PEGFP-C1


Clontech

pEGFP-C1 *NOBAMURA B2*

Product Quick Links:


- [Annotated Map/MCS \(pdf\)](#)
- [Sequence \(text\)](#)
- [Sequence \(pdf\)](#)
- [Restriction Digest \(pdf\)](#)
- [Description](#)
- [Use](#)
- [Feature Locations](#)
- [Primer Sites](#)
- [Propagation](#)
- [References](#)
- [License Statement](#)

- PRODUCTS**
 - eCatalog
 - Product Families
 - Technologies
 - New Products
 - Sales Reps
 - Special Offers
- NEWS**
 - Clontechniques
 - Highlights
 - Coming Soon
 - Conferences
- SUPPORT**
 - User Manuals
 - Vector Information
 - Citations
 - Application Notes
 - FAQs
 - Email Tech Support
- LITERATURE**
 - Printed Literature
 - Brochures
 - Clontechniques
 - Order Forms
 - Product Lists
- ABOUT**
 - BD Biosciences
 - Contact Us
 - Employment



[Click here to view the MCS for this vector.](#)

Restriction Map and Multiple Cloning Site of pEGFP-C1. (Unique restriction sites are in bold.) The *Xba* I and *Bcl* I sites (*) are methylated in the DNA provided by CLONTECH. If you wish to digest the vector with these enzymes, you will need to transform the vector into a dam- host and make fresh DNA.



MADE by: CHRIS SULLIVAN
8/1/03
Subcloned KpnI-APAF
fragment from pGUA3.1.PURO
NOB B2 INTO SAME
SITES.

- CONFIRMED BY
ZFA.



Description

pEGFP-C1 encodes a red-shifted variant of wild-type GFP (1-3) which has been optimized for brighter fluorescence and higher expression in mammalian cells. (Excitation maximum = 488 nm; emission maximum = 507 nm.) pEGFP-C1 encodes the GFPmut1 variant (4) which contains the double-amino-acid substitution of Phe-64 to Leu and Ser-65 to Thr. The coding sequence of the EGFP gene contains more than 190 silent base changes which correspond to human codon-usage preferences (5). Sequences flanking EGFP have been converted to a Kozak consensus translation initiation site (6) to further increase the translation efficiency in eukaryotic cells. The MCS in pEGFP-C1 is between the EGFP coding sequences and the SV40 poly A. Genes cloned into the MCS will be expressed as fusions to the C-terminus of EGFP if they are in the same reading frame as EGFP and there are no intervening stop codons. SV40 polyadenylation signals downstream of the EGFP gene direct proper processing of the 3' end of the EGFP mRNA. The vector backbone also contains an SV40 origin for replication in mammalian cells expressing the SV40 T-antigen. A neomycin resistance cassette (*neo^r*), consisting of the SV40 early promoter, the neomycin/kanamycin resistance gene of Tn5, and polyadenylation signals from the Herpes simplex thymidine kinase (HSV TK) gene, allows stably transfected eukaryotic cells to be selected using G418. A bacterial promoter upstream of this cassette expresses kanamycin resistance in *E. coli*. The pEGFP-C1 backbone also provides a pUC origin of replication for propagation in *E. coli* and an f1 origin for single-stranded DNA production.

>gb|AF174533.1| **G** Nodamura virus RNA1 segment, complete sequence
Length=3204

Score = 740 bits (820), Expect = 0.0
Identities = 413/415 (99%), Gaps = 0/415 (0%)
Strand=Plus/Plus

```
Query 85 ATGACAAACATGTCATCGGCTTACGAGCTAATCAAGTCACATTCAGCCAAAGCTGGAGCAG 144
|||||
Sbjct 2744 ATGACAAACATGTCATCGGCTTACGAGCTAATCAAGTCACATTCAGCCAAAGCTGGAGCAG 2803

Query 145 CTGGCTCAGGAGACGCAAGCAACGATCCAAACGGTCAATGATCGGGATCCCCAACGTCAAC 204
|||||
Sbjct 2804 CTGGCTCAGGAGACGCAAGCAACGATCCAAACGGTCAATGATCGGGATCCCCAACGTCAAC 2863

Query 205 AAGGATCTGCGAGCGTTCTGCGAGTTCCGTGACCCGTACAGCACCCAGCGGGCGTATCGAGCG 264
|||||
Sbjct 2864 AAGGATCTGCGAGCGTTCTGCGAGTTCCGTGACCCGTGACCCAGCACCCAGCGGGCGTATCGAGCG 2923

Query 265 ACGAACAGCCCTGCTCATCAAACCGGAGTCCGAGCAGCGCTTCGCGGGAGGAGCTGGAC 324
|||||
Sbjct 2924 ACGAACAGCCCTGCTCATCAAACCGGAGTCCGAGCAGCGCTTCGCGGGAGGAGCTGGAC 2983

Query 325 CTGGCGAGCGGACGTCGCCGCCCGGTCGCCAGCTAAACAACAGCTGGCGGAGCTC 384
|||||
Sbjct 2984 CTGGCGAGCGGACGTCGCCGCCCGGTCGCCAGCTAAACAACAGCTGGCGGAGCTC 3043

Query 385 GAGATGGAAATCAAGCCAGGGCACCAACAGTGGCCCAAGTAAGCCGCAGCCGGAAGGCC 444
|||||
Sbjct 3044 GAGATGGAAATCAAGCCAGGGCACCAACAGTGGCCCAAGTAAGCCGCAGCCGGAAGGCC 3103

Query 445 GCAGCCGACGTCCTCCGTGGCCAGCTGGGTGCGGTGGCCGTGTAATGAGTGAT 499
|||||
Sbjct 3104 GCAGCCGACGTCCTCCGTGGCCAGCTGGGTGCGGTGGCCGTGTAATGAGTGAT 3158
```

arrows indicate 2 silent point mutations