

# pDG\* One-Step Cloning Protocol

**!** Oligos for guide insertion into plasmid need to be of the following form:  
 5' -CACCGNNNNNNNNNNNNNNNNNNNNNN-3'  
 3' -CNNNNNNNNNNNNNNNNNNNNNNCAA-5'

And:

5' -ACCGNNNNNNNNNNNNNNNNNNNNNGT-3'  
 3' -NNNNNNNNNNNNNNNNNNNNNNCAAAT-5'

**!** If the first N on the top strand for each guide is a G, it should be excluded.

- ?** The N's in the top strands comprise the guide sequences, which target the matching gRNA binding sequences followed by PAMs in the genomic DNA.
- ?** The overhangs allow the oligos to match the cut sites in the pDG\* plasmids caused by *BbsI* digestion.
- ?** The extra G/C in the first pair of oligos completes the U6 promoter, this doubles as the first base in the guide which is why it should be excluded if the guide starts with a G.
- ?** The extra GT/CA in the second pair of oligos completes the gRNA scaffold.

1. Mix the following reagents in a **PCR tube** for each of the two inserts:

Reagent	Amount
MQ H <sub>2</sub> O	6.5 µL
NEB T4 DNA Ligase Reaction Buffer (10x)	1 µL
bottom oligo (100 µM)	1 µL
top oligo (100 µM)	1 µL
NEB T4 Polynucleotide Kinase (10 U/µL)	0.5 µL
<b>Total</b>	<b>10 µL</b>

**?** T4 Polynucleotide Kinase phosphorylates the 5' end of each oligo to allow more efficient ligation when cloning into a plasmid.

**!** 2. Place each of the two mixtures in a thermocycler with the following parameters:

<b>1</b>	37 °C	30 min
<b>2</b>	95 °C	5 min
<b>3</b>	Ramp to 25 °C @ 0.1 °C/s	∞

**?** 37 °C for optimal T4 PNK activity, 95 °C to melt then ramp down to anneal.

3. Dilute the 2 sets of phospho-annealed oligos 1:125 with **MQ H<sub>2</sub>O** in a **1.5 mL tube**.

Reagent	Amount
MQ H <sub>2</sub> O	124 µL
phospho-annealed oligo	1 µL
<b>Total</b>	<b>125 µL</b>

4. Mix the following reagents in a **PCR tube**:

Reagent	Amount
MQ H <sub>2</sub> O	11.5 µL
pDG* plasmid (100 ng/µL)	1 µL
phospho-annealed oligo 1 (1:125)	1 µL
phospho-annealed oligo 2 (1:125)	1 µL
NEBuffer 2.1 (10x)	2 µL
DTT (10 mM)	1 µL
ATP (10 mM)	1 µL
NEB <i>BbsI</i> (5 U/µL)	1 µL
NEB T4 DNA Ligase (400 U/µL)	0.5 µL
<b>Total</b>	<b>20 µL</b>

- ❓ Dithiothreitol (DTT) is a reducing agent used to prevent the formation of disulfide bonds in the Ligase, allowing it to maintain its function.
- ❓ ATP is integral to the function of the reaction that the Ligase catalyses.

🕒 5. Place in thermocycler with the following parameters:

1	37 °C	5 min
2	16 °C	5 min
3	Go to step 1	5 times
4	4 °C	∞

- ❓ Cycle between optimal temperature for *BbsI* which cuts the plasmid open and T4 DNA Ligase (16 °C) which ligates the phospho-annealed oligos into those cut sites and destroys the *BbsI* sites. Alternatively, the original pDG\* plasmid is re-assembled by the Ligase without insertion of the phospho-annealed oligos, in which case the plasmid is cut again in the next 37 °C step. 5 cycles is sufficient to get enough plasmid containing the custom phospho-annealed oligos.

- ➡ Transform [Invitrogen Subcloning Efficiency DH5α Competent Cells](#) with 5 μL of mixture, following the recommended protocol.
- ➡ Verify insertion of both guides are present in the plasmid by one of the following:
  - Digesting with any restriction enzymes that have new binding sites generated by the insertion of the guide.
  - Digesting with *BbsI* and another restriction enzyme present in the backbone.
  - Sequencing both inserts with **seq gRNA U6 V1**: 5'-GGTTTCGCCACCTCTGACTTG-3' and **bgh PA F**: 5'-TGCATCGCATTGTCTGAGTAGG-3'.