

Protocol for cloning MHC ligand minigenes into PresentER vector

Ron Gejman, Scheinberg Lab (ron.gejman@gmail.com). Updated 4/13/2018

Reagents:

- Phusion High Fidelity Polymerase (NEB)
- SfiI restriction enzyme 20,000 units/ml (NEB)
- Calf Intestinal Phosphatase 10,000 units/ml (NEB)
- T4 DNA Ligase 400,000 units/ml (NEB)
- Primers:
 - T7_SfiI: CGACTCACTATAGGGCCGTATTGGCC
 - T3_SfiI: AGTGATTTCCGGCCTGTTTGGCC
- Qiagen MinElute Kit (to PCR purify amplicons < 100nt)
- Gel purification kit (Qiagen or other)
- PresentER Cassette (GFP or mCherry version)

Oligonucleotide design:

Order standard desalted oligonucleotide(s) with the following design:
GGCCGTATTGGCCCGCCACCTGTGAGCGGG + Peptide DNA + TAAGGCCAAACAGGCC

For example:

RMFPNAPYL:

GGCCGTATTGGCCCGCCACCTGTGAGCGGGAGGATGTTTCCTAACGCGCCCTA
CCTGTAAGGCCAAACAGGCC

ALYVDSLFFL:

GGCCGTATTGGCCCGCCACCTGTGAGCGGGGCTCTCTATGTGGACTCTTTATT
TTTCCTTAAGGCCAAACAGGCC

1. Resuspend oligonucleotides at 100µM stock.
2. Make 1:10,000 dilutions of the oligo (~1.5 ng/ul):
3. Amplify oligos in 50µl reactions:
 - 10 µl Phusion HF Buffer
 - 1 µl dNTPs (freshly thawed from an aliquot)
 - 2.5 µl of 10µM T7_SfiI primer
 - 2.5 µl of 10µM T3_SfiI primer
 - 0.5 µl Phusion DNA polymerase
 - 32.5 µl water
 - 1 µl oligo DNA

33 cycles w/ 72°C annealing and run 5µl out on a 2.2% gel to check for amplification.
(10s denaturation; 30s annealing; 30s extension; 5 minute final extension)

Amplicon is ~100nt

4. PCR purify with MinElute kit. Elute in 10-20µl EB.
5. Digest amplicon and PresentER Cassette with SfiI at 50°C for 2h - o/n
- 6.

Amplicon restriction digest

8µl amplicon DNA
2µl cutsmart
1µl SfiI
9µl water

20µl total reaction volume

Check amplicon digest on a 2.2% gel. Digest should result in a faint band at ~70nt (for a 8-9aa MHC-I ligand) and elimination or reduction in the ~100nt amplicon band. It helps to use ethidium bromide gels because it is more sensitive than SYBR safe.

Purify digested amplicon with MinElute kit and elute in 10-20µl of water. Technically some digestion products will be smaller than the lower limit of the MinElute kit, but in practice this hasn't been a problem.

PresentER Cassette (Vector) restriction digest

5µg PresentER Cassette Plasmid DNA
5µl Cutsmart
1µl SfiI
X µl water

50µl total reaction volume

Digest 20-50µg of DNA total in reactions of 5µg (50µl) each.

Check digest on a gel to verify it has gone to completion. If digest is incomplete, add 1µl SfiI + 1µl cutsmart + 8µl water to each reaction and digest for another hour.

Add 1µl of Calf intestinal phosphatase to each reaction and incubate 1h at 37°C

Gel purify digested vector.

6. Ligate digested+purified amplicon and digested+CIP treated+gel purified vector at 16°C overnight
 - x µl MinEluted insert
 - X µl backbone
 - 2µl 10x T4 buffer
 - 1µl T4 ligase (400 units)
 - X µl water
 -
 - 20µl total volume

5 moles Insert : 1 mole Vector

9. Transform 10µl of NEB STABLE cells with 2µl of ligation product, pick colonies; miniprep;
10. Sanger sequence with custom primer CATCAGTCCGGGTACCTAC or universal primer GAG (GTTCGACCCCGCCTCGATCC)