


# OVERVIEW


## Step 1: Landing site targeting



### Summary

- Electroporate sgRNA/Cas9 expression plasmid and single-stranded oligodeoxynucleotide (ssODN)
- Dilute to clonal density after 72 hours
- Isolate colonies and screen by PCR


## Step 2: Array plasmid integration



### Summary

- Electroporate integrase expression plasmid and array targeting vector
- Drug selection or FACS sorting
- Dilute to clonal density
- Isolate colonies and screen by PCR

## Optional Step 3: Selection cassette removal



### Summary

- Electroporate Cre or Flp expression plasmid
- Dilute to clonal density after 72 hours
- Isolate colonies and screen by PCR

## attP ssODN Templates

$(n)_{70}$  <sup>XbaI</sup> tctagacc<sup>PhiC31 attP</sup>ccaactggggtaaccttgagttctctcagttggggg<sup>XhoI</sup>ctcgag $(n)_{70}$   
 5' homology Bxb1 attP 3' homology  
 $(n)_{65}$  tctagagtcgtggtttgtctggtcaaccaccg<sup>Bxb1 attP</sup>cggtctcagtggtgtacggtacaaccccgactcgag $(n)_{65}$

## Recombination Arms

### PhiC31 attL

ggtgccaggggcgtgcccttgagttctctcagttggggg

### PhiC31 attR

ccaactggggtaaccttgggctccccgggcgcg

### Bxb1 attL


tcggccggcctgtgcagcagcggcgggtctcagtggtgtacggtacaaccccgac

### Bxb1 attR


gtcgtggtttgtctggtcaaccaccg<sup>Bxb1 attR</sup>cggtctccgctcgtcaggatcatccgggc

## PCR Genotyping

### Step 1



### Step 2



## Protocols for Insertion of cuO-/tetO-arrays into Mouse ESCs

### Reagents

0.1% gelatin -- Sigma ES-006-B  
 Accutase Solution -- Thermo Fisher MT25058CI  
 0.05% trypsin solution -- Thermo Fisher 25300-054  
 ES quality FBS -- Thermo Fisher SH3007002E  
 Neon Electroporation Kit -- Thermo Fisher MPK1096  
 GoTaq PCR Mix -- Promega M7123  
 D-PBS w/o Ca<sup>++</sup> and Mg<sup>++</sup> -- Thermo Fisher 14190-144  
 DMEM/F12 -- Thermo Fisher 11320-033  
 Neurobasal -- Thermo Fisher 21103-049  
 N2 Supplement -- Thermo Fisher 17502-048  
 B27 w/ Retinoic Acid -- Thermo Fisher 17504-044  
 7.5% BSA -- Thermo Fisher 15260-037  
 100X GlutaMax -- Thermo Fisher 35050-061  
 1-thioglycerol – Sigma M6145  
 PD03259010 – Selleckchem 1036  
 CHIR99021 – Selleckchem S2924  
 LIF – Peprotech 250-02

### Recipes

**2i media** -- 50% DMEM/F12:50% Neurobasal, N2 Supplement, B27 Supplement, 0.05% BSA, 2 mM GlutaMax, 150 uM 1-thioglycerol, 1 uM PD03259010, 3 uM CHIR99021, 10<sup>6</sup> U/L LIF

**Washing Media** -- DMEM, 0.05% BSA

### Oligonucleotide attP Homology Template Design

Design and order an oligonucleotide from IDT with the PhiC31 or Bxb1 attP sequence, flanking restriction sites for genotyping (XbaI and XhoI are recommended), and 65-70 bases of homology to the targeted genomic region of interest on each side of the attP insert (for a total of 200 bp). This can be ordered from IDT as an Ultramer. Place the attP insert sequence such that it is located between the CRISPR protospacer recognized by the sgRNA and the protospacer adjacent motif (PAM). This will ensure the knock-in allele will not be re-edited by CRISPR/Cas9. For more details regarding designing CRISPR sgRNAs and generating CRISPR/sgRNA expression vectors from pX330, see <http://www.addgene.org/crispr/zhang/>.



### Landing Site Targeting

#### Electroporation

- 1\_Culture mouse ES cells for at least one passage subsequent to a thaw.
  - 2\_Add 0.5 mL 0.1% gelatin to 24 well plate for desired number of electroporation conditions. Incubate for 30 minutes.
  - 3\_Remove 0.1% gelatin and add 0.5 mL 2i media without antibiotics. Place in 37°C tissue culture incubator.
  - 4\_Dissociate ES cells with Accutase Solution for roughly 2 minutes at 37°C.
  - 5\_Transfer cell suspension to 15 mL conical and bring volume up to 10 mL with Washing Media.
  - 6\_Centrifuge at 300xg for 3 minutes at room temperature.
  - 7\_Wash cell pellet once with 10 mL Washing Media.
  - 8\_Centrifuge at 300xg for 3 minutes at room temperature.
  - 9\_Resuspend cells at approximately 10<sup>6</sup> cells/mL in Washing Media. Perform cell count.
  - 10\_Transfer 1.5 x 10<sup>5</sup> cells per electroporation condition to a fresh 15 mL conical.
- \_\_\_ # conditions x 150k cells = \_\_\_ total # of cells needed

- 11\_Bring up to 10 mL in Washing Media.
  - 12\_Centrifuge at 300xg for 3 minutes at room temperature.
  - 13\_Wash once with 10 mL D-PBS w/o Ca<sup>++</sup> and Mg<sup>++</sup>.
  - 14\_Centrifuge at 300xg for 3 minutes at room temperature.
  - 15\_Carefully remove all traces of supernatant.
  - 16\_Resuspend cell pellet in R solution (see Neon Electroporation kit) at 10<sup>7</sup> cells/mL.  
This works out to be 15 uL R solution for each 150,000 cells (i.e. each well of a 24-well plate). At this point, the protocol should be streamlined and completed as quickly as possible.
  - 17\_Aliquot 15 uL of cell suspension into individual 1.5 mL tubes, one for each electroporation condition (i.e. each well of a 24-well plate).
  - 18\_Add DNA for electroporation to each 1.5 mL tube.  
400 ng pX330-sgRNA plasmid (at ~500-1000 ng/uL)  
1 uL of 10 uM attP ssODN homology template
  - 19\_Use 10 uL Neon pipette tip to electroporate cells according the manufacturer's protocol. After electroporation, immediately add cells to one well of 24-well plate with pre-warmed 2i media.  
Neon pipette will aspirate 10 uL (100,000 cells) for electroporation. Each condition is prepared at 15 uL volume to ensure that no air bubbles are introduced into Neon pipette, which will inhibit electroporation.
- Recommended electroporation conditions**  
1400 V, 10 ms pulse width, 3 pulses
- 20\_Allow 72 hours for recovery after electroporation, then proceed to dilution to clonal density (listed below).

## Array Plasmid Integration

### *Electroporation*

- 1\_Culture mouse ES cells harboring attP knock-in for at least one passage subsequent to a thaw.

### **Steps 2-17 (as listed above)**

- 18\_Add DNA for electroporation to each 1.5 mL tube.

#### **PhiC31- or Bxb1-mediated Plasmid Integration**

400 ng integrase expression plasmid (at ~500-1000 ng/uL)  
400 ng attB donor plasmid with tetO/cuO array (at ~500-1000 ng/uL)

#### **PhiC31/Bxb1-mediated Double Plasmid Integration**

300 ng PhiC31 integrase expression plasmid (at ~500-1000 ng/uL)  
300 ng Bxb1 integrase expression plasmid (at ~500-1000 ng/uL)  
300 ng PhiC31 attB donor plasmid with tetO array (at ~500-1000 ng/uL)  
300 ng Bxb1 attB donor plasmid with cuO array (at ~500-1000 ng/uL)

- 19\_Use 10 uL Neon pipette tip to electroporate cells according the manufacturer's protocol. After electroporation, immediately add cells to one well of 24-well plate with pre-warmed 2i media.  
Neon pipette will aspirate 10 uL (100,000 cells) for electroporation. Each condition is prepared at 15 uL volume to ensure that no air bubbles are introduced into Neon pipette, which will inhibit electroporation.

#### **Recommended electroporation conditions**

1400 V, 10 ms pulse width, 3 pulses

- 20\_Allow 72 hours for recovery after electroporation.
- 21\_Passage electroporated cells. If using drug selection, add selective drug(s) to 2i media without antibiotics.  
Appropriate concentration of drug should be determined empirically for each batch. However, a sensible place to start is listed below.  
G418 selection – 500 ug/mL final concentration  
Puromycin – 1 ug/mL final concentration

- 22\_If using drug selection plasmids, culture cells under selection for the time necessary to kill all non-expressing cells. If using fluorescent reporter plasmids, culture cells for 7 days to allow unintegrated plasmids to be lost.  
Necessary time for selection should also be determined empirically by treating non-expressing ES cells. However, the durations we observed were the following:  
G418 selection – 6-7 days  
Puromycin – 2 days

- 23\_After selection is complete, proceed to dilution to clonal density (listed below).

## Selection Cassette Removal

### *Electroporation*

- 1\_Culture mouse ES cells with integrated cuO and/or tetO arrays for at least one passage subsequent to a thaw.

### **Steps 2-17 (as listed above)**

- 18\_Add DNA for electroporation to each 1.5 mL tube.

#### **FRT cassette removal**

400 ng Flp expression plasmid (at ~500-1000 ng/uL)

**Floxed cassette removal**

400 ng Cre expression plasmid (at ~500-1000 ng/uL)

**Double removal of FRT and Floxed cassette**

400 ng Flp expression plasmid (at ~500-1000 ng/uL)

400 ng Cre expression plasmid (at ~500-1000 ng/uL)

19\_ Use 10 uL Neon pipette tip to electroporate cells according the manufacturer's protocol. After electroporation, immediately add cells to one well of 24-well plate with pre-warmed 2i media.

Neon pipette will aspirate 10 uL (100,000 cells) for electroporation. Each condition is prepared at 15 uL volume to ensure that no air bubbles are introduced into Neon pipette, which will inhibit electroporation.

**Recommended electroporation conditions**

1400 V, 10 ms pulse width, 3 pulses

**Dilution to Clonal Density**

1\_ Add 7 mL 0.1% gelatin to 10 cm dish for each condition to isolate clones. Incubate for 30 minutes.

2\_ Dissociate ES cells with Accutase Solution for roughly 2 minutes at 37°C.

3\_ Transfer cell suspension to 15 mL conical and bring volume up to 10 mL with Washing Media.

4\_ Centrifuge at 300xg for 3 minutes at room temperature.

6\_ Resuspend in 100 uL Washing Media. Perform cell count.

7\_ Seed the pre-gelatinized 10 cm dish with 8 mL 2i media with ~5000 cells.

8\_ Allow single cells to grow into large, round colonies, changing media every 2 days.

Single ES cells in 2i will take 5-7 days to grow into large colonies ideal for picking.

**Isolation of ES Cell Colonies**

1\_ Treat a 96-well plate with 50 uL 0.1% gelatin per well. Incubate for 30 minutes.

2\_ Remove gelatin. Add 25 uL 0.05% trypsin to each well.

3\_ On day of colony picking, transfer 10 cm dish to tissue culture microscope for colony picking. Wipe down microscope and place in tissue culture hood if possible.

4\_ While viewing under microscope, use pipette tip to gently dislodge attached colony and aspirate colony into pipet tip. Transfer to one 96-well with trypsin.

Many colonies may already be dislodged at this point and can be directly aspirated from the growth media.

5\_ Continue collecting colonies and fill the 96-well plate column-wise (adding clones to wells A1-H1, followed by A2-H2). Once the second column is complete, add 100 uL 2i + 5% ES-quality serum to first column to quench trypsinization. Pipet up and down to dissociate colonies.

6\_ Repeat step 5 for each column sequentially (i.e. add 2i + 5% ES-quality serum to second column after third column colonies are picked).

17\_ After all colonies are picked, place in incubator and grow until near confluent (~2-3 days).

Media should be changed to 100 uL 2i on the following day to remove any trace trypsin activity.

18\_ Split plate of clones into 2 replica 96-well plates. Use one replicate for extracting DNA to use for PCR screening. Maintain the second plate in culture or freeze until genotyping is complete.

For additional details on these steps, we recommend *Manipulating the Mouse Embryo: A Laboratory Manual*

## PCR Screening

### *attP Knock-in*

1\_Design genomic PCR primers that flank the attP knock-in site and have no overlap with the oligonucleotide template designed above.

2\_Set up PCR reaction using GoTaq PCR Mix

#### **PCR Reaction (per sample)**

H <sub>2</sub> O	8.5 uL
2x GoTaq Mix	10 uL
10 uM Primer Mix	0.5 uL
Genomic DNA (10-25 ng/uL)	1 uL
-----	
Total	20 uL

- 1- 95°C -- 5 minutes
  - 2- 95°C -- 30 seconds
  - 3- Primer Tm -- 30 seconds
  - 4- 72°C -- 45-60 seconds
- Repeat steps 2-4 an additional 34 times
- 5 - 72°C -- 5 minutes
  - 6 - Hold at 12°C

3\_Run 10 uL product on agarose gel. Identify samples that have an upward shifted band that is consistent with attP insertion.

4\_To the remaining 10 uL of PCR product from promising samples, add 10 units XhoI or XbaI restriction enzyme.  
Ensure that the enzyme used is not expected to digest the PCR product from unmodified genomic DNA.

5\_Incubate at 37°C for 2 hours.

6\_Run digested PCR on agarose gel. Identify samples where upper band has been digested and verify restriction fragment size. These samples are attP insertion clones.

These clones can be further validated by sequencing.

### *Array Plasmid Integration*

1\_Use the same PCR primers used to identify attP knock-in clones.

2\_Set up PCR reaction using GoTaq

#### **PCR Reaction (per sample)**

H <sub>2</sub> O	6.125 uL
2x GoTaq Mix	7.5 uL
10 uM Primer Mix	0.375 uL
Genomic DNA (10-25 ng/uL)	1 uL
-----	
Total	15 uL

- 1- 95°C -- 5 minutes
  - 2- 95°C -- 30 seconds
  - 3- Primer Tm -- 30 seconds
  - 4- 72°C -- 45-60 seconds
- Repeat steps 2-4 an additional 34 times
- 5 - 72°C -- 5 minutes
  - 6 - Hold at 12°C

3\_Run products on agarose gel. Identify samples where the upward shifted band from the parental attP knock-in line has disappeared due to failure to amplify large cuO or tetO array insertion.

These clones can be further validated by using primer pairs that span the attR and attL junctions to confirm correct integration.

## References

Nagy, A., Gerstenstein, M., Vintersten, K., Behringer, R., (2003). *Manipulating the Mouse Embryo: A Laboratory Manual* (Cold Spring Harbor Laboratory Press).