

# Short-homology-mediated CRISPR/Cas9-based method in fission yeast

Hayashi & Tanaka 2019, G3 <https://doi.org/10.1534/g3.118.200976>

## Two plasmids (Cas9 + gRNA)

Transform w/ Cas9 vector  
 ↓ + thiamine, 4 days  
 Transform w/ gRNA vector + donor  
 ↓ - thiamine, 6 days  
 Streak transformants on YES plate  
 ↓ 2~3 days  
 Confirm the genome editing  
 e.g. phenotype  
 ↓  
 colony PCR  
 ↓  
 sequencing

## Single plasmid (gRNA-Cas9)

Transform w/ gRNA-Cas9 vector + donor  
 ↓ - thiamine, 6 days  
 Streak transformants on selective plate  
 ↓ - thiamine, 2 days  
 Streak cells on YES plate  
 ↓ 2~3 days  
 Confirm the genome editing  
 e.g. phenotype  
 ↓  
 colony PCR  
 ↓  
 sequencing

## 1. Design gRNA oligomers

### 1-1. Design gRNAs

Choose gRNA sequence near the position you would like to introduce mutation/insertion. Or, go to the following web sites (e.g. CRISPRdirect) to find the appropriate gRNA candidates. You can input ~500 bp genomic sequence, select the genome of organisms (e.g. fission yeast) and initiate a search. Avoid gRNA that has low GC content (less than 40%).

### 1-2. Check the gRNA specificity

Choose gRNA that has low number of similar sequences in the genome. Avoid using gRNA that has many A/T sequences near PAM.

### 1-3. order oligomers for gRNA cloning

gRNA - Fw 5' -caccXXXXXXXXXXXXXXXXXXXX-3'  
 gRNA - Rv 5' -aaacXXXXXXXXXXXXXXXXXXXX-3'

*Do not include PAM (NGG)!*

## Web tool for gRNA design

CRISPRdirect

<https://crispr.dbcls.jp>

E-CRISP

<http://www.e-crisp.org/E-CRISP/designcrispr.html>

Bähler lab genome Regulation

<http://bahlerweb.cs.ucl.ac.uk/cgi-bin/crispr4p/webapp.py>

## 2. Cloning gRNA

### 2-1. annealing of oligos

CACCXXXXXXXXXXXXXXXXXXXX  
XXXXXXXXXXXXXXXXXXXXCAA


*10 x annealing buffer	3 ul
gRNA Fw (100 uM)	6 ul
gRNA Rv (100 uM)	6 ul
dH <sub>2</sub> O	15 ul /30 ul

Mix the oligos in a PCR tube. Heat at 95°C 2 min in thermal cycler then leave at R.T for 30 min~1h.

Or you can use a program for annealing (95°C 2 min, cool-down to 25°C (- 2°C/ min) of thermal cycler.

\*10 x annealing buffer: 0.1 M Tris pH7.5~8, 0.5 M NaCl, 10 mM EDTA

### 2-2. digestion gRNA expression plasmid w/ *BbsI*



After plasmid DNA preparation, heat for 30 min at 65°C and keep in a freezer.  
Cas9 plasmid is less stable, it is desirable to keep at -80°C.

Digestion

Incubation at 37°C, 1 h

10 x buffer	10 ul
Plasmid (1 ~ug)	X ul
<i>BbsI</i> (NEB) (10 U/ ul)	0.4 ul
dH <sub>2</sub> O	X ul /100 ul

Check linearization of plasmid by gel electrophoresis

Heat-inactivation of *BbsI* (65°C, 20 min)

EtOH precipitation

(If works, you can purify using column DNA purification kit. Sometimes you might get a low recovery because of large-sized plasmid.)

Keep at -80°C.

### 2-3. ligation of annealed oligo and gRNA expression plasmid

e.g. (TOYOBO Ligation high ver.2)

Incubation at 16°C, 30 min	Ligation high (1/2 volume of DNA)	1 ul
	<i>BbsI</i> digested Plasmid (~10 ng/ul)	1 ul
	1/100 diluted annealed oligo	1 ul / 3 ul

## 2-4. transformation

Add ligated plasmid (~5 ul) to 50 ul of competent cells (e.g. *DH5α*)  
Incubate cells on ice, for 30 min  
Heat at 42°C, 30 sec  
Cool down on ice, for 2 min  
Spread cells on LB-Amp plate  
Incubate overnight at 37°C


## 2.5. gRNA insertion check

Grow single colony in LB-Amp broth at 37°C overnight and isolate plasmid DNA next day using column DNA preparation kit. Heat plasmid DNA at 65°C for 30 min and keep in a freezer.

For check gRNA insertion, do colony PCR with gRNA Fw and the proper Rv primer if you like.  
And/or gRNA insertion is confirmed by sequencing.

## 3. Design donor oligomers

### 3. 1. for introduced point mutation



Find the cleavage site that is located at - 3 bp from PAM. Choose 40~50 bp of sequence near the cleavage site (20~25 bp of homologous sequences from the cleavage site). The knock-in frequency is very high if the mutation site is located near the cleavage site.

Order donor oligomer (~50 bp). If you would like to introduce mutation at a high frequency, use mixed (Fw and Rv oligomer) or double stranded oligomer for transformation.

As donor oligos, using single stranded oligo (Fw) is able to introduce genome editing. Mixed and/or double stranded oligo are precisely introduced the knock-in at high frequency.


### 3. 2. for knock-in tag

Mixed or double stranded oligo are better for the introduction of precise insertion.

### 3. 2. for knock-in tag (continued)



Order oligomers (~25 bp homology + first ATG + tag + ~25 bp homology).



## 4. Transformation w/ CRISPR/Cas9 plasmid and donor oligo

### - two plasmids protocol -

You can stock the yeast cells carrying Cas9 expression plasmid at -80°C freezer. For stock, the cells carrying Cas9 expression plasmid are cultured in EMM – uracil + thiamine medium to repress Cas9 expression. Before transformation, the Cas9 expressing cells are cultured in EMM – uracil – thiamine medium for overnight.

### - single plasmid protocol -

The yeast cells are cultured in EMM5 – thiamine medium for overnight.

Transformation: e.g. *S.pombe* transformation kit (Wako) (Li-Ac method)

Log phase culture ( $1\sim 2 \times 10^7$  cells/ml) is harvested and suspended with the medium to adjust cell number to  $1\sim 2 \times 10^9$  cells/ml.

e.g. harvest 10 ml of  $1\sim 2 \times 10^7$  cells/ml culture and suspend w/ 100~200 ul of medium.

Mixed:	~1 ug of gRNA expression plasmid (or gRNA-Cas9 expression plasmid)
	donor DNA (single stranded oligo: 1 nmol, double stranded oligos: 900 pmol)
	carrier DNA ( e.g. 2 ul of 2 mg/ml salmon testes DNA)
	$1\sim 2 \times 10^7$ cells (10 ul)
	Transformation reagent (including PEG)(45 ul)

Incubate the mixed cells at 37°C for 2~3 h. If host cell is temperature-sensitive, you can incubate the cells at permissive temperature (e.g. for 6h~overnight at 25°C).

Spread the cells gently on selective medium plate

Incubate plate at 32°C for 6~7 days. If the cell would be obtained temperature sensitive phenotype by genome editing, incubate at permissive temperature (e.g. 10~ days at 25°C).

## 5. Confirmation of mutation/insertion introduction

### 5.1. choose small transformants

#### - two plasmids protocol -

Choose small transformants and streak transformants on YES plate to get single colonies.  
Incubate plate for 2~3 days.

*do not choose big colonies!*

#### - single plasmid protocol -

Choose small transformants and streak transformants on selective medium plate and incubate plate for 2~3 days. Then, streak the cells on YES plate to get single colonies.

### 5.2. phenotype check

If your mutant has some phenotypes, you can select the mutants before you confirm the genome editing by sequencing.

### 5.3. plasmid loss check

Streak single colonies on both of YES and selective medium plate to check the loss of plasmids from the cells. If you use *ura4<sup>+</sup>* -marked plasmid, you can use 5-FOA plate to remove the plasmids.

### 5.4. direct colony PCR

Choose colonies that lost the plasmids for direct colony PCR.

e.g.	10 x PCR Buffer		1 ul
	2 mM dNTPs		0.8 ul
	primers (2.5 uM)	Fw & Rv	0.4 ul each
	Taq polymerase (5 U/ul)		0.05 ul
	dH <sub>2</sub> O		X ul
	*cells		/10 ul

\*touch fresh colony with 200 ul tip and mix well w/ reaction mixture in a PCR tube

94°C 2 min - 94°C 30 sec. - 52 °C 40 sec.- 72°C (1 kb/ 1 min) x 36 cycles - 72°C 10 min  
- hold 12 °C

### 5.5. Electrophoresis of PCR product

Electrophoresis of PCR products (1 ul) to check the size/ rearrangement/ insertion.

### 5.6. sequencing

Remove primers and dNTPs from PCR product. You can purify PCR product DNA using purification kit. Or, use ExoSAP (ThermoFisher).

e.g.	PCR product	2.5 ul	
	ExoSAP-IT	0.25 ul	
	<i>express</i>		
	dH <sub>2</sub> O	0.75 ul	/1 sample

37°C 4 min - 80°C 1 min - hold 4°C

Use 1 ul of treated PCR product for sequencing analysis.