

Dilute this reaction 1:1000. This is your insert.

Step 3: Digest the shuttle vector plasmid:

x μ L	6 μ g of DNA
6 μ L	10x Cutsmart Buffer
2 μ L	Bsal-HF enzyme
(52-x) μ L	Water
60 μ L	Total volume

Digest at 37°C overnight (16-18 hours at least). Gel extract to obtain your shuttle vector. Do not CIP-treat. Ensure the vector concentration after elution is less than 50ng/ μ L (dilute if necessary).

Step 4: Ligate annealed oligos into the Shuttle vector:

1 μ L	Linearized Shuttle vector
2.5 μ L	gRNA annealed primers (1:1000 dilution from above)
2 μ L	NEB 10x T4 DNA Ligase Buffer
1 μ L	T4 DNA Ligase
13.5 μ L	Water
20 μ L	Total Volume

Incubate at room temperature for 20 minutes, then transform 25 μ L competent cells (*our lab uses homegrown competent cells*) with 5 μ L of ligation reaction (*This can be changed to maintain a 1:5 ratio of ligation mix to competent cells*). Add 37.5 μ L SOC medium (*1.5x volume of competent cells used*) and recover for 15 minutes at 37°C. Plate full volume on LB+Amp plates. Incubate plates at 37°C overnight.

(*For first use after digesting the Shuttle vector, perform a “Shuttle vector only” negative control to test for digestion efficiency. Make sure you don’t get any colonies. Repeat digestion if necessary.*)

Step 5: Amplification of HH-gRNAf+e-HDV precursor:

Pick individual colonies and swirl in 10 μ L of SOC (in a PCR strip tube). Set up the following reactions:

Reaction 1

2.5 μ L	Swirled Culture
5 μ L	EcoRI-HH-<Gene>-Fwd
5 μ L	HDV-NotI-Rev
2.5 μ L	Taq Buffer
1 μ L	dNTPs
1.5 μ L	DMSO
0.5 μ L	Taq Polymerase
32 μ L	Water
50 μ L	Total Volume

In parallel, set up the following colony PCR to verify successful integration of the protospacer:

Reaction 2

1 μ L	Swirled Culture
2 μ L	gRNA-TOP (10 μ M dilution)
2 μ L	T7 seq Reverse
1 μ L	Taq Buffer
0.4 μ L	dNTPs
0.6 μ L	DMSO
0.2 μ L	Taq Polymerase
12.8 μ L	Water
20 μ L	Total Volume

Run 2 μ L of Reaction 1 and 20 μ L of Reaction 2 on a gel. Reaction 1 should yield a ~215bp product, while Reaction 2 should yield a 690bp product (if the protospacer was successfully integrated into the shuttle vector). PCR purify the remaining 48 μ L of Reaction 1 (for colonies that have a positive band for Reaction 2) and elute in 52 μ L. Set up the following digestion reaction:

52 μ L	PCR-purified product
6 μ L	10x Cutsmart Buffer
1 μ L	EcoRI-HF enzyme
1 μ L	NotI-HF enzyme
60 μ L	Total volume

Digest for 2 hours at 37°C. Gel extract. This is your “EcoRI-HH-gRNA-HDV-NotI” insert.

Step 6: Clone HH-gRNAf+e-HDV precursor into single CRISPR plasmid:

Set up the following digestion reaction:

x μ L	6 μ g of single CRISPR plasmid
6 μ L	10x Cutsmart Buffer
1 μ L	EcoRI-HF enzyme
1 μ L	NotI-HF enzyme
(52-x) μ L	Water
60 μ L	Total volume

Digest at 37°C overnight. Gel extract to obtain the CRISPR destination vector. Set up ligation at 1:3 molar ratio using this vector and the insert obtained in Step 5.

x μ L	Linearized CRISPR destination vector
y μ L	Digested and purified HH-gRNA-HDV insert (from Step 5)
2 μ L	NEB 10x T4 DNA Ligase Buffer
1 μ L	T4 DNA Ligase
(20-x-y) μ L	Water
20 μ L	Total Volume

Ligate for 2 hours at room temperature. Transform 25 μ L competent cells (*our lab uses homegrown competent cells*) with 5 μ L of ligation reaction (*This can be changed to maintain a 1:5 ratio of ligation mix to competent cells*). Add 37.5 μ L SOC medium (*1.5x volume of competent cells used*) and recover for 1 hour at 37°C. Plate and incubate at 37°C overnight.

The next day, validate colonies by colony PCR using Citrine-ColonyPCR-Fwd and 10 μ M dilution of gRNA-BOT primers. Validate insert by sequencing and electroporate chick embryos at 2.5 μ g/ μ L concentration.