

# A rapid sgRNA cloning protocol for Multicolor CRISPR (V0315)

This protocol was adapted from Ma H *et al*, PNAS, 2015. Hanhui Ma (Thoru Pedeson's lab @ UMMS).

## A. Design sgRNA oligos

Fwd (all): ACCGNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN

Rev (Sp): AAACNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN

Rev (Nm): CAACNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN

Rev (St): AGACNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN

PAM Choice:

Sp dCas9: NGG

Nm dCas9: NNNNGATT

St dCas9: NNAGAAW

## B. Annealing oligos

### Reagents:

	Stock	ml
• Annealing Buffer		
100 mM Tris-HCl pH8.0	1 M	0.25
50 mM NaCl	5 M	0.25
1 mM EDTA	0.5 M	0.05
<u>ddH<sub>2</sub>O</u>		<u>24.5</u>
Total		25

- 100 μM oligos

Annealing:	μl
Annealing buffer	40
Fwd oligo (100 μM)	5
<u>Rev oligo (100 μM)</u>	<u>5</u>
Total	50

- Incubate @ 95 °C for 3 min and slowly cool down to room temperature;
- Dilute 5 μl of annealed oligos to 245 μl water and final concentration is 200 nM (about 3ng/ul, ready for use).

## C. Cloning annealed oligos into vectors

### Reagents:

- pLH-sgRNA vectors\*
- Enzymes and Buffer in reaction
- Stbl3 competent cells
- LB- Ampicillin plates

Reaction mix:	μl
pLH-sgRNA vector (100 ng/μl)	1
10XSmartCut Buffer (NEB)	1
10 mM ATP (NEB)	1
BbsI (NEB)	0.5
T7 DNA ligase (NEB)	0.3
ddH <sub>2</sub> O	5.2
<u>Annealed oligos (200 nM)</u>	<u>1</u>
	10

- Incubate @ 37 °C for 15 min;
- Transform 5 μl of reaction mix into Stbl3 competent cells and spread on LB-Amp plates.

## D. Minipreps and Sequencing

### Reagents:

- QIAprep Spin Miniprep Kit
- Sequencing primer pLKO.1-Rseq:  
CTATTCTTTCCCCTGCACTGTACCC

### \* Notes of pLH-sgRNA vectors:

pLH-sgRNA2 for Sp dCas9.

pLH-nmsgRNA1.1 for Nm dCas9.

pLH-st1sgRNA1.1, pLH-st1sgRNA2.1 and pLH-st1sgRNA3.1 for St1 dCas9.

All these vector plasmids are available @

[Addgene](#) and they require CcdB Survival cells

for growth.