


I. Cloning of sgRNAs into CRISPR/Cas9-T2A-Reporter vectors

pX330-T2A-Reporter (mCherry/BFP)



5' – **CACCG**NNNNNNNNNNNNNNNNNNNNNNNNNN – 3'
 3' – **C**NNNNNNNNNNNNNNNNNNNNNNNNNN**CAAA** – 5'

1. Digest pX330-T2A-mCherry/BFP

Reagents	Volume
pX330-T2A-Reporter	1 μ g
Fastdigest BbsI (fermentas)	1 μ l
FastAP (Fermentas)	1 μ l
10x FastDigest buffer	2 μ l
H ₂ O	Up to 20 μ l
Digest for 60 min at 37°C 70°C for 5min	

2. Gel purify digested pX330-T2A-mCherry/BFP using QIAquick gel extraction kit and elute in H₂O

3. Phosphorylate and anneal each pair of oligos

Reagents	Volume
Oligo 1 (stock 100 μ M)	1 μ l (10 μ M)
Oligo 2 (stock 100 μ M)	1 μ l (10 μ M)
10x T4 ligation buffer (NEB)	1 μ l
H ₂ O	6.5 μ l
T4 PNK (NEB)	0.5 μ l
Total	10 μl

Anneal parameters	
37°C	30 min
95°C	5 min and then ramp down to 25°C at 0.1°C/second, hold 4°C

Annealed oligo duplex (10 μ M) will be diluted 1:20 to get 500nM stock analysed with 4% agarose gel (used 5 μ l). From 500nM stock will be further diluted 1:10 get 50nM oligo duplex for cloning. Oligo duplex are stored at -20°C

4. Ligation reaction

Reagents	Volume
BbsI digested pX330-T2A-mC/BFP	Final 50ng
Phosphorylated and annealed duplex (1:200)	1.5 μ l


dilution) (step 3)	
10X T4 ligation buffer (NEB)	1 μ l
T4 Ligase (NEB, M)	0.5 μ l
H2O	Up to 10 μ l
Incubate at 22°C for 15 min	

5. Transformation

6. Plasmid isolation

7. Sequencing to confirm the insert

II. Cloning of new inserts into sgRNA-pX330-T2A-Reporter vectors



1. Digest pX330-sgRNA-T2A-mCherry/BFP

Reagents	Volume
pX330-sgRNA-T2A-Reporter	1 μ g
Fastdigest NheI (fermentas)	1 μ l
Fastdigest EcoRI (fermentas)	1 μ l
FastAP (Fermentas)	1 μ l
10x FastDigest buffer	2 μ l
H2O	Up to 20 μ l
Digest for 60min at 37°C 70°C for 5min	

2. Gel purify digested pX330-sgRNA-T2A-mCherry/BFP using QIAquick gel extraction kit and elute in H₂O

3. Digest pX330-T2A-Reporter-P2A-Ad4 E1B55K/E4orf6 vector

Reagents	Volume
The vector	1 μ g
Fastdigest NheI (fermentas)	1 μ l
Fastdigest EcoRI (fermentas)	1 μ l
10x FastDigest buffer	2 μ l
H2O	Up to 20 μ l
Digest for 60min at 37°C 70°C for 5min	

4. Gel purify digested Reporter-P2A-Ad4 E1B55K/Ad4 E4orf6 fragments using QIAquick gel extraction kit and elute in H₂O

5. Ligation reaction


Reagents	Volume
NheI/EcoRI digested vector	Up to 50ng
NheI/EcoRI digested fragment	Up to 30ng
10X T4 ligation buffer (NEB)	1 μ l
T4 Ligase (NEB, M)	0.5 μ l
H2O	Up to 10 μ l
Incubate at 22°C for 15 min	

6. Transformation

7. Plasmid isolation

8. Sequencing to confirm the insert

III. Cloning of shRNAs into CRISPR/Cas9-T2A-reporter-hH1-shRNA vectors



shRNA oligos

Sense

Loop

Antisense

5' -gatccNNNNNNNNNNNNNNNNNNNNNNNTCAAGAGNNNNNNNNNNNNNNNNNNNNNNNtttttc -3'
3' - gNNNNNNNNNNNNNNNNNNNNNNNAGTTCTCNNNNNNNNNNNNNNNNNNNNNNNaaaagaatt-5'

1. Digest pX330-T2A-mCherry/BFP-hH1 vector

Reagents	Volume
pX330-T2A-mCherry/BFP-hH1 vector	1µg
Fastdigest BamH1 (fermentas)	1µl
Fastdigest AflIII (fermentas)	1µl
FastAP (Fermentas)	1µl
10x FastDigest buffer	2µl
H2O	Up to 20µl
Digest for 60min at 37°C 70°C for 5min	

2. Gel purify digested pX330-T2A-mCherry/BFP-hH1 using QIAquick gel extraction kit and elute in H₂O

3. Phosphorylate and anneal each pair of oligos

Reagents	Volume
Oligo 1 (stock 100µM)	1µl (10µM)
Oligo 2 (stock 100µM)	1µl (10µM)
10x T4 ligation buffer (NEB)	1µl
H2O	6.5µl
T4 PNK (NEB)	0.5µl
Total	10µl

Anneal parameters	
37°C	30 min
95°C	5 min and then ramp down to 25°C at 0.1°C/second, hold 4°C

Annealed oligo duplex (10µM) will be diluted 1:20 to get 500nM stock analysed with 4% agarose gel (used 5µl). From 500nM stock will be further diluted 1:10 get 50nM oligo duplex for cloning. Oligo duplex are stored at -20°C

4. Ligation reaction

Reagents	Volume
pX330-T2A-mCherry/BFP-hH1 vector	Final 50ng
Phosphorylated and annealed duplex (1:200 dilution) (step 3)	1.5 μ l
10X T4 ligation buffer (NEB)	1 μ l
T4 Ligase (NEB, M)	0.5 μ l
H2O	Up to 10 μ l
Incubate at 22°C for 15 min	

5. Transformation

6. Plasmid isolation

7. Sequencing to confirm the insert