

For cloning the spacer (target site) sequence onto the 5' end of single guide RNAs (sgRNAs) for human and bacterial expression vectors:

1) Human sgRNA expression vectors (U6 promoter)




Plasmid	Addgene ID	Description	Components	Cut with	spacer
BPK1520 & MLM3636*	65777 & 43860	<i>S.pyogenes</i> sgRNA cloning	U6-BsmBIIcassette-Sp.sgRNA	BsmBI	17-20
BPK2301	65778	<i>S.thermophilus1</i> sgRNA cloning	U6-BsmBIIcassette-St1.sgRNA	BsmBI	20-23
BPK2660 & VVT1**	70709 & 65779	<i>S.aureus</i> sgRNA cloning	U6-BsmBIIcassette-Sa.sgRNA	BsmBI	21-23

Ordering oligonucleotides to generate the spacer (target) sequence for human sgRNA vectors:

For BPK1520, BPK2301, BPK2660, & VVT1:


*For MLM3636:



* The N₂₀ sequence can be of variable length, as indicated in the table above

** VVT1 contains the full crRNA/tracrRNA sequence, whereas BPK2660 is a truncated and more active sgRNA

2) Bacterial sgRNA/Cas9 expression vectors (T7 promoter)




Plasmid	Addgene ID	Description	Components	Cut vector with
BPK764	65767	<i>S.pyogenes</i> Cas9/sgRNA cloning vector	T7-SpCas9-T7-Bsal-Sp.sgRNA	Bsal
MSP712*	65768	<i>S.pyogenes</i> <u>d</u> Cas9/sgRNA cloning vector	T7-Sp <u>d</u> Cas9-T7-Bsal-Sp.sgRNA	Bsal
MSP1673	65769	<i>S.thermophilus1</i> Cas9/sgRNA cloning vector	T7-St1Cas9-T7-BspMI-St1.sgRNA	BspMI
BPK2101	65770	<i>S.aureus</i> Cas9/sgRNA cloning vector	T7-SaCas9-T7-Bsal-Sa.sgRNA	Bsal

Ordering oligonucleotides to generate the spacer (target) sequence for bacterial sgRNA vectors:

For BPK764, MSP1673, and BPK2101:

*For MSP712:



Cloning steps (see following page for protocol):

- Order oligos** with indicated overhangs (*notes: 1. N₂₀ corresponds to your target of interest for the top strand oligo, use the reverse complement N₂₀ for the bottom strand oligo. 2. The N₂₀ sequence can be extended or truncated, should different spacer lengths be desired. 3. Please note the overhang differences for MLM3636 & MSP712 vs. the other vectors*)
- Anneal and phosphorylate (optional)** the top & bottom strand oligos for a given spacer pair
- Digest and dephosphorylate (optional)** the vector backbone, into which annealed oligos will be cloned (*note: different IIS enzymes depending on vector of choice*)
- Ligate** the phosphorylated oligo duplexes into the corresponding cut backbone
- Transform** the ligation product into *E.coli* (*note: human sgRNA vector cloning should be plated on ampicillin or carbenicillin, while bacterial sgRNA cloning should be plated on chloramphenicol*)
- Prep and sequence** the resulting clones

Cloning protocol:

1) **Suspend oligos:** to 100 μM in 0.1xTE buffer, dilute 1:10 to 10 μM in H_2O

2) **Generate phosphorylated (optional) oligo duplex:**

1 μL top strand oligo (10 μM)
1 μL bottom strand oligo (10 μM)
2 μL 10x T4 DNA Ligase buffer (NEB)
0.5 μL T4 PNK (NEB) (optional)
16 or 15.5 μL H_2O
20 μL reaction volume

Reaction conditions

- Incubate at 37°C for 60 minutes (only if phosphorylating)
- Heat to 95°C for 5 minutes
- Cool to 10°C at -5°C/minute

3) **Digest & dephosphorylate (optional) the vector backbone:**

x μL plasmid (2-3 μg plasmid DNA)
4 μL 10x buffer (corresponding to enzyme)
y μL restriction enzyme
0.5 μL CIP (NEB) (optional)
z μL H_2O (to sum to 40)
40 μL reaction volume

Reaction conditions

- 37°C or 55°C for 4-16 hours

<u>Enzyme</u>	<u>Buffer</u>	<u>Temp</u>
BsmBI	NEB 3.1	55°C
Bsal-HF	Cutsmart	37°C
BspMI	NEB 3.1	37°C

- Isolate cut plasmid

4) **Ligation:**

x μL digested plasmid from step 3 (10-40 ng)
1 μL of oligo duplex from step 2 (0.5 μM)
1 μL 10x T4 DNA ligase buffer (NEB)
0.5 μL T4 DNA ligase (NEB)
y μL H_2O
10 μL reaction volume

Reaction conditions

- Incubate at 24°C for 10 minutes
- Chill on ice prior to transformation

5) **Transform the ligation into *E.coli* and plate on appropriate selective media:**

- **ampicillin** or **carbenicillin** for human sgRNA expression vectors
 - **note:** *high copy number vectors*
- **chloramphenicol** for bacterial Cas9/sgRNA expression vectors
 - **note:** *low copy number vectors* – appropriate liquid media recommended when starting O/N cultures

6) **Sequence prepped plasmids:**

- For human sgRNA vectors, use a primer 5' of the U6 promoter in the pUC19 backbone
 - For example, the primer αS280 : 5'-CAGGGTTATTGTCTCATGAGCGG-3'
- For bacterial Cas9/sgRNA expression vectors, use a reverse primer 3' of the sgRNA
 - For example, the T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'