

Materials:

Plasmids/DNA:	CROPseq-iT7 (Addgene Plasmid #211699) Addgene Brunello pooled library #73178
Oligonucleotides:	
cropseq_fwd	TATCTTGTGGAAAGGACGAAACA
cropseq_rev	GCCTTATTTTAACTTGCTATTCTAG
MiSeq_gRNA_fwd_S0	ACACTCTTTCCCTACACGACGCTCTTCCGATCTATCTTGTGGAAAGGACGAAACAC
MiSeq_gRNA_fwd_S1	ACACTCTTTCCCTACACGACGCTCTTCCGATCTCATCTTGTGGAAAGGACGAAACAC
MiSeq_gRNA_fwd_S2	ACACTCTTTCCCTACACGACGCTCTTCCGATCTGCATCTTGTGGAAAGGACGAAACAC
MiSeq_gRNA_fwd_S3	ACACTCTTTCCCTACACGACGCTCTTCCGATCTAGCATCTTGTGGAAAGGACGAAACAC
MiSeq_gRNA_fwd_S4	ACACTCTTTCCCTACACGACGCTCTTCCGATCTCAACATCTTGTGGAAAGGACGAAACAC
MiSeq_gRNA_fwd_S6	ACACTCTTTCCCTACACGACGCTCTTCCGATCTTGACCATCTTGTGGAAAGGACGAAACAC
MiSeq_gRNA_fwd_S7	ACACTCTTTCCCTACACGACGCTCTTCCGATCTACGCAACATCTTGTGGAAAGGACGAAACAC
MiSeq_gRNA_fwd_S8	ACACTCTTTCCCTACACGACGCTCTTCCGATCTGAAGACCCATCTTGTGGAAAGGACGAAACAC
MiSeq_gRNA_rev	TGACTGGAGTTCAGACGTGTGCTCTTCCGATCTCCGACTCGGTGCCACTTTTTTCAA
Barcoding primer fwd	AATGATACGGCGACCACCGAGATCTACAC X ₁₀ ACACTCTTTCCCTACACGACGCT
Barcoding primer rev	CAAGCAGAAGACGGCATACGAGAT X ₁₀ GTGACTGGAGTTCAGACGTGTGCT
Enzymes/Buffers:	FastDigest Esp3I (ThermoFisher, Cat.No. FD0454) FastDigest Clal (ThermoFisher, Cat.No. FD0143) 10x FastDigest Buffer (ThermoFisher, Cat.No. B64) FastAP (ThermoFisher, Cat.No. EF0651) DTT NEBNext® High-Fidelity 2X PCR Master Mix (New England Biolabs, Cat.No. M0541L) 2x Gibson Assembly MasterMix (New England Biolabs, Cat.No. E2611L)
Bacteria/Media:	Endura™ Competent Cells DUOs (BiosearchTechnologies Cat.No. 60242-1) Endura Recovery Medium LB medium LB agar plates Ampicillin (final concentration 100 µg/mL)
Materials/Kits:	QIAprep 2.0 Spin Miniprep Columns (Cat.No. 27115) Qiagen buffer PE (Cat. No. 19065) Qiagen buffer PB (Cat. No. 19066) Zymo DNA Clean & Concentrator® -25 (Cat.No. D4005) BioRad cuvette 0.1 cm (Cat.No. #1652089) PureLink™ HiPure Plasmid Filter Maxiprep Kit (ThermoFisher, Cat.No. 210017)

Protocol:**Library cloning**

- Digest vector plasmid CROPseq-iT7 with Esp3I and Clal
 - Prepare a master mix for three reactions containing:
 - 9 µL 10X FastDigest buffer
 - 3 µg backbone plasmid
 - 4.5 µL Esp3I FastDigest
 - 4.5 µL Clal FastDigest
 - 4.5 µL Fast AP
 - 0.9 µL 100 mM DTT
 - Ad 90 µL with H₂O
 - Prepare 30 µL reactions in PCR stripes and put in thermocycler at 37°C for 2.5 hours
 - Take 2 µL aliquots for running on an agarose gel, pool reactions, purify using a Qiagen PCR cleanup column (blue column), elute in 30 µL H₂O, and quantify using a NanoDrop
 - store at 4°C
- PCR Brunello sgRNA library
 - Prepare a master mix for 16 reactions containing:
 - 150 µL 2x NEBNext HighFidelity PCR MasterMix
 - 12 µL Lenti-guide Puro Brunello sgRNA plasmid library (10 ng/µL)
 - 30 µL mix of primers cropseq_fwd and cropseq_rev, 5 µM each
 - 108 µL H₂O
 - Prepare 25 µL reactions in PCR stripes and run the following program in a thermocycler:
 - 98°C 30 sec
 - 98°C 10 sec
 - 55°C 10 sec
 - 72°C 15 sec
 - 72°C 2 min
 - 12°C hold
 - Take 2 µL aliquots for running on an agarose gel, pool reactions after incubation, purify using a QIAprep 2.0 Spin Miniprep Column, pool eluate (3x 25 µL), run 2 µL of purified DNA on an agarose gel, and quantify pooled DNA using a NanoDrop
 - Store at 4°C

- Gibson Assembly
 - Set up a master mix for 10 reactions:
 - 19 μ L digested backbone
 - 10 μ L amplified sgRNA insert
 - 100 μ L 2X Gibson Assembly master mix
 - 71 μ L H₂O
 - Prepare 20 μ L reactions and incubate at 50°C for 1h
 - Pool 10 reactions and concentrate using the Zymo DNA Clean & Concentrator kit
 - Elute in 25 μ L H₂O, measure concentration using a NanoDrop

Amplification of library - Electroporation of Endura competent cells

- Prepare eight 0.1 cm cuvettes (BioRad) and place the cuvettes on ice to cool down
- Set the electroporator to exponential mode with 1800 V, 10 μ F, 600 Ohms
- Thaw four vials (50 μ L each) of Endura competent cells (BioSearch Technologies) on ice (~10-20 minutes)
- Per cuvette, prepare a 1.5 mL Eppendorf tube with 975 μ L Endura recovery medium at room temperature
- Once the cells are thawed, prepare 25 μ L aliquots in a PCR stripe on ice and add 2 μ L (60.2 ng) library, mix by stirring in tube with the pipet tip – do not pipet up and down, since this may warm up the cells
- Transfer the mix to the gap of the pre-cooled cuvettes, tap cuvette on the bench a few times to ensure distribution over the whole gap, dry metal contacts, and quickly electroporate before water is condensing on the contacts again
- If electroporation was successful, transfer bacteria to one of the prepared tubes with recovery medium
- Allow recovery period of all electroporated tubes for 1 hour at 37°C, 400 rpm
- In the meantime, pre-warm 6 LB-Agar plates and prepare a 1L flask with 500 mL LB-Amp
- After recovery period, pool library samples in separate tubes and mix by inverting
- Prepare dilutions for plating, library A-D
 - A: 10 μ L plus 990 μ L LB (final: 1:100)
 - B: from (A) 10 μ L plus 90 μ L LB (final: 1:1,000)
 - C: from (A) 10 μ L plus 990 μ L LB (final: 1:10,000)
 - D: from (B) 10 μ L plus 990 μ L LB (final: 1:100,000)
 - plate 80 μ L of A-D on LB Agar-Amp plates and incubate at 37°C o/n
- Pool remaining library dilution (A) with pooled culture, transfer to 1L flask with medium, and incubate in big shaker at 37°C, 230 rpm o/n

Calculation of library diversity:

- Count the colonies on the agar plates with less than 500 colonies, e.g.:
Plate B (1:1,000) – 220 colonies
Plate C (1:10,000) – 21 colonies
Plate D (1:100,000) – 2 colonies
- Multiply with the count of each plate with the dilution factor from the corresponding plate:
 - B: 220 colonies x 1,000 = 220,000
 - C: 21 colonies x 10,000 = 210,000
 - D: 2 colonies x 100,000 = 200,000
- Multiply with the relative amount of dilution plated from the start culture (80 μ L of 8 mL \rightarrow 100x) to estimate the total number of clones:
 - B: 100 x 220,000 = 22,000,000 total clones
 - C: 100 x 210,000 = 21,000,000 total clones
 - D: 100 x 200,000 = 20,000,000 total clones

\rightarrow Estimated library complexity after electroporation between 20 and 22 Mio \rightarrow mean of all three plates calculated ~21 Mio.

\rightarrow For a library with 80.000 sgRNAs each guide is present ~263 times

- Purify the DNA from maxiprep culture with PureLink™ HiPure Plasmid Filter Maxiprep Kit according to manufacturer's instructions and dissolve the DNA pellet in 300-1000 μ L H₂O

Validation via Sequencing

- Prepare a 10 ng/μL dilution of the library maxiprep
- Prepare a primer mix consisting of the eight staggered forward (fwd) primers MiSeq_sgRNA_S0-S8 and the reverse (rev) primer MiSeq_sgRNA_rev
 - Each fwd primer should have a final concentration of 0.625 μM, the rev primer should have a final concentration of 5 μM
 - Mix 2 μL of each fwd primer with 16 μL of the rev primer
 - Label as "5 μM Brunello staggered primer mix"
- Prepare a master mix for PCR1 for 4 replicates, each at a total volume of 20 μL:
 - 10 μL 2x NEBNext HighFidelity PCR MasterMix
 - 2 μL 5 μM Brunello staggered primer mix
 - 1 μL DNA
 - 7 μL H₂O
- Prepare reactions in PCR stripes and run the following program in a thermocycler:
 - 98°C 3 min
 - 98°C 20 sec
 - 65°C 20 sec
 - 72°C 30 sec
 - 72°C 3 min
 - 4°C hold

| x20
- Prepare a master mix for PCR2, each at a total volume of 20 μL:
 - 10 μL 2x NEBNext HighFidelity PCR MasterMix
 - 2 μL unique fwd / rev barcoding primer mix (5 μM each)
 - 2 μL PCR1
 - 6 μL H₂O
- Run thermocycler using the same program as for PCR1
- Pool reactions after incubation, purify using a QIAprep 2.0 Spin Miniprep Column, pool eluate, run 2 μL of purified DNA on an agarose gel, and quantify pooled DNA using a NanoDrop
- Run a quality control of the sequencing library with ~10,000 reads on an Illumina MiSeq
- If QC looks good, run the library with at least 100,000,000 reads on an Illumina NextSeq 2000

Analyse Library Coverage:

- In web browser (Firefox or Google Chrome) open <https://www.jsb-lab.bio/LibCounter.htm>
- In the blank field "Library" paste the reference library copied from a text file consisting of two columns separated by tab:
 - name of the targeted gene
 - sequence of the sgRNA
- Load the sequencing FASTQ files below at "choose files"
- Press "Start" to run the analysis and retrieve a count table of sgRNAs present in your library
- Sort the list of sgRNAs by reads per sgRNA from lowest to highest read count
- Calculate the fold change between the 10th and 90th percentile as a measure of even coverage
- A fold change below 15 is considered sufficient to perform genome-scale screens
- If the library quality is sufficient, proceed with production of lentivirus

Lentivirus production

Day 0: Seed HEK cells

Seed 15x10⁶ HEK293T cells in DMEM 10% FCS + Ciprofloxacin per 15 cm dish

Day 1: transfection

- Prepare Lipofectamine 2000 in OptiMEM, per 15 cm dish mix:
 - 120 μL Lipofectamine 2000 + 1.66 mL OptiMEM
 - Incubate 5 minutes at room temperature
- Prepare a DNA mix in OptiMEM, per 15 cm dish mix:
 - 7,040 ng pMD2.G
 - 10,560 ng psPAX2
 - 14,080 ng CROPseq-iT7-Brunello plasmid library
 - Fill up with OptiMEM to 1.78 mL
- Combine both mixes and incubate for 20 minutes at room temperature
- Add the mix dropwise to the HEK cells (go in a spiral pattern from the outside of the dish to the center)
- After 4-6 hours, change the medium in each dish to 30 mL fresh DMEM supplemented with 10% FCS and 10 μg/mL Ciprofloxacin

Day 3: harvest

- After 48 hours incubation, harvest the lentivirus
- Pool supernatant from all dishes in 50 mL falcon tubes
- Clear lentiviral supernatant by filtering through a 0.45 μm filter using a syringe
- Aliquot in Sarstedt screw cap 1.5 mL tubes and flash freeze in liquid N₂
- Store in labeled box at -80°C