

## Production of retroviruses using FuGENE-6

3/21/00

- 1) On the day before transfection, plate  $0.5-1 \times 10^6$  293T cells onto 6cm dishes in 4-5 ml DME/10%IFS. The cells are ready for transfection after 18-20 hours, or when they are about 25-30% confluent. Plate slightly more cells when making VSV-G pseudotyped viruses. For 10 cm dishes, double the quantities given here.
- 2) Pipette 6 $\mu$ l FuGENE into serum-free medium (at room temperature) so that the final volume after the addition of the DNA will be 100 $\mu$ l. DO NOT TOUCH THE TUBE WITH UNDILUTED FUGENE.
- 3) Add 2 $\mu$ g total DNA (in TE). For instance, if using pCL vectors (Ampho, 10A1, or Eco; Naviaux (1996) *J. Virol*, vol 70, p 5701), use 1 $\mu$ g transfer vector (pBabe, pMSCV-IRES-GFP, etc.) and 1  $\mu$ g of the pCL vector. For making VSV-G retroviruses by triple transfection use 0.9 $\mu$ g gag/pol expression vector, 0.1  $\mu$ g VSV-G expression vector, and 1  $\mu$ g transfer vector.
- 4) Do not vortex; swirl to mix.
- 5) Incubate the FuGENE/medium/DNA solution at room temperature for 15 minutes.
- 6) Drip the mixture onto the 293T cells using a Pipetman.
- 7) Incubate overnight at 37 C.
- 8) On the day following transfection feed the cells with 3-4 ml fresh medium.
- 9) Harvest viral supernatants at 48 and 72 hours (or 36 and 60 hours). Ecotropic and GaLV viral titers at 72 hours will be approximately 20-40% less than at 48 hours. VSV-G viral titers are essentially the same at 48, 72, and 96 hours. Filter with a 45  $\mu$ M syringe filter.
- 10) In a typical viral infection protocol for adherent cells: use 3.5 ml of 36-48 hour viral sup, with 8  $\mu$ g/ml polybrene, to infect  $1-2 \times 10^6$  cells (20-30% confluent) on a 10 cm dish. After 3 hours remove the sup and add fresh medium.
- 11) Repeat the infection 24 hours later with the 60-72 hour sup. 24 hours after the second infection, split the cells; apply drug selection or sort the day following the split. Although two infections are optimal, one infection is often enough.
- 12) Please occasionally take the time to sing happy songs about viruses.

### Notes about using FuGENE to make retroviruses:

- Do not let FuGENE come into contact with any plastic surface without diluting first in serum free medium (do not aliquot; store in original polypropylene tubes).
- Antibiotics do not seem to impair the transfection or viral titers.
- Whether FuGENE or DNA is added first to the medium doesn't seem to matter. It's OK to add the two or three plasmids separately.
- For collecting VSV-G viruses, use medium buffered with 10-25 mM HEPES to limit pH-dependent cell killing by the VSV-G protein.
- With this protocol we have made ECO viruses with titers of  $10^7$ /ml and VSV-G, 10A1, and AMPHO viruses with titers of  $10^6$ /ml.
- FuGENE is made by Roche, formerly Boehringer Mannheim. Cat # 1 814 443, 1 ml, \$275, about \$1.65 per transfection.

Scott Dessain and Hua-yin Yu  
Weinberg Lab