

Protocol for Cloning of single guide RNAs for Cas9 genome editing using LentiCRISPR

Materials

- Tango Buffer
- Oligos (standard synthesis, purification)
- Screw cap tubes or Eppendorf tubes (e.g. 1.5 ml)
- CRISPR target plasmid with BsmBI insert site (eg LentiCRISPR-EGFP, mCherry etc...)
- T4 DNA ligase (NEB, M0202M or M202L, **make aliquots of T4 DNA ligase buffer and use only freshly thawed aliquots for reaction**)
- Esp3I restriction enzyme (Thermo Scientific, ER0451; **NOT NEB BsmBI**, BsmBI requires 55°C for efficient cleavage – not compatible with the reaction here)
- Forward test primer for colony PCR: 5'-GACTATCATATGCTTACCGT-3' (specific for U6 promoter in LentiCRISPR)
- Mini- and/or Midi-Prep kits (e.g. Qiagen)

Preparation

- screen your target genomic locus for CRISPR guide RNA binding sites (e.g. <http://crispr.mit.edu>)
- use the target designer excel sheet to get the forward and reverse primers for your sgRNA
- paste the 20 nucleotide sequence of your favourite guide RNA(s) (**with PAM**) into the target designer generator
- order the generated forward and reverse oligos (lyophilized or at a concentration of 100 µM, desalted is fine)

Day 1 – Annealing of single-stranded oligos and Ligating of double-stranded oligos into LentiCRISPR

- for lyophilized oligos elute in an appropriate volume of H₂O to obtain an oligo concentration of 100 µM
- pipette the following reaction for each pair of forward and reverse oligos into a screw cap tube and close tightly:

10 µl oligo forward
10 µl oligo reverse
10 µl Tango Buffer (or similar)
70 µl H₂O

- bring about up to 95°C for 5 mins
- Let oligos return to room temperature over 1-2h

Inserting annealed oligos into LentiCRISPR

- set up following reaction for ligation of annealed double-stranded oligos into LentiCRISPR:

150 ng of LentiCRISPR plasmid
1 µl of double-stranded oligos
2µl NEB T4 ligase buffer
to 18 µl with H₂O
0.5 µl Esp3I
0.5 µl T4 ligase

- run the following reaction in a thermocycler

37°C 5 min
16°C 10 min 10 cycles (can be extended to 50 cycles over night)
37°C 15 min
80°C 5 min

- Note: There is an additional BsmBI(Esp3I) site in the tagBFP sequence in the LC-BFP vector. For this vector the final 37°C step should be removed from the cycle protocol to improve yield

University Children's Hospital Zurich 2016
Leukemia Lab, Bornhauser Cell Death Group

- transform ligation mix into competent bacteria (electroporation or heat shock)
- plate transformed bacteria on LB-Agar + Ampicillin

Day 2 – Identify bacterial clones with LentiCRISPR + double-stranded oligo insert

- run a colony PCR for 3-5 colonies per transformed LentiCRISPR (3 colonies is almost always enough as this procedure should have 80-100% efficiency)
- Colony PCR program:


95°C 3 min
95°C 30 sec
55°C 30 sec
72°C 10 sec 25 cycles
72°C 10 min

- prepare separate LB-agar Amp plate with a number for each clone that will be picked
- prepare PCR reaction mix (20 µl) with the primers

Forward: Forward Test Primer (GACTATCATATGCTTACCGT)

Reverse: reverse single-stranded guideRNA oligo that was used for annealing of double-stranded oligos (day 1)

- pick single colonies with toothpick or pipette tip and streak on separate numbered plate and subsequently submerge toothpick/pipette tip in one well of PCR reaction buffer
- after picking all colonies, remove toothpicks/pipette tips, run PCR and put numbered plate with streaks at 37°C for 2-5h
- load PCR reactions on 2% agarose gel
 - o clones with LentiCRISPR + double-stranded oligo insert should give a band around 100 bp
 - o clones with empty LentiCRISPR should give no band



- pick positive clone(s) from numbered plates to start overnight cultures for mini or midi preps
- You can also sequence plasmid using the forward vector if you wish to confirm insert
- Proceed to Lentiviral production using ps.Pax2 and p.VSV.G in 293T Cells
- Some hints for Lentiviral work:
 - o LentiCRISPR should only be handled by those with working knowledge and safety training in the production and use of Lentivirus.
 - o LentiCRISPR is known to have very low viral yields due to the size of the plasmid.
 - o We have improved our viral production by increasing the amount of viral plasmid transfected into 293T cells
 - o We use 4ug LentiCRISPR + 1.5 ug ps.Pax2 + 1.5ug p.VSV.G in 3mL plated in a 6 well dish
 - o We have had success using both PEI and CaCl₂/HEBS transfections
 - o We change the media at 6 hours after transduction, then collect viral supernatant at 24 and 48 hours.
 - o Don't expect to get 50-100% transduced cells using raw supernatant. If you wish to get high transduction you will need to concentrate.
 - o We generally aim for 5% transduction of cells and then sort from there.