

Ple236 Promoter PCR (pEMS1208)**Pleiades Promoter Project**

MiniPromoter: Ple236
pEMS#: 1208
Expected product size (bp): 377

Reaction components	Vol/Rxn (µl)
H ₂ O	15.15
10X PCR buffer*	2.5
50 mM MgCl ₂ *	0.75
2.5 mM dNTPs**	2
10 µM forward primer	1.25
10 µM reverse primer	1.25
Taq Pol. (5 U/µl)*	0.1
DNA***	2
Total Volume of Rxn:	25

* Taq Polymerase set from Invitrogen (Cat no.18038-042)

** dNTPs from Invitrogen (Cat no.10297-018)

***Approximately 100 ng DNA used for samples, approximately 5 ng used for plasmid control

Samples run on a 2% agarose gel (containing SYBRsafe (Invitrogen Cat no. S33102))

Cycling conditions:	Step	Temp	Time	Note
	1	94°C	3 min	
	2	94°C	1 min	
	3	61°C	1 min	
	4	72°C	45 sec	repeat steps 2-4 34 times
	5	72°C	5 min	
	6	4°C	hold	

Primers:

Name	Sequence (5' to 3')	T_m (°C)	Notes
oEMS2926	AACCTGGGTCACGTGGACTACC	56.5	Sense primer for Ple236, in region "3"
oEMS2885	GCTTTGTCCCCAGACTGGA	56.2	Antisense primer for Ple236, in region "Prom"