


Plasmid DNA Description:

The synthetic gene SPIKE_EUK_18S was assembled from synthetic oligonucleotides and/or PCR products. The fragment was inserted into pMA-T. The plasmid DNA was purified from transformed bacteria and concentration determined by UV spectroscopy. The final construct was verified by sequencing. The sequence congruence within the insertion sites was 100%. See the accompanying data sheets for sequences and find the original ABI trace files as well as the assembled sequences electronically on disk. 5 µg of the plasmid preparation were lyophilized for shipping.

Plasmid Map:



invitrogen
by Thermo Fisher Scientific

Quality Assurance Documentation: 16ABQ5JP

Ref. No.: 1917357

Designation: E. coli K12 OmniMAX™ 2 T1R

Gene name: SPIKE_EUK_18S

Gene size: 212 bp

Vector backbone: pMA-T

Cloning sites: SfiI / SfiI

Quantity: ~5 µg Plasmid DNA

Note: Please dissolve lyophilized DNA in 50 µl distilled water or 10 mM Tris-HCl (pH 8.0). We recommend sequence verification after each transformation step.

Date: 13 June 2016

Dr. Norbert Greifenegger

Quality control

www.thermofisher.com

geneartsupport@thermofisher.com