

## Open Source Wnt Nomenclature

- Wnt1-STOP: Wnt1 sequence with a STOP codon cloned into the Entry Vector, pENTR/D-TOPO
- Wnt1-NONSTOP: Wnt1 sequence without a STOP codon cloned into the Entry vector pENTR/D-TOPO
- Wnt1 (untagged): Wnt1 sequence with a STOP codon cloned into the Destination vector, pcDNA3.2/V5-DEST. These Wnts do NOT express the V5 tag found in the backbone.
- Wnt1-V5 (tagged): Wnt1 sequence without a STOP codon cloned into the Destination vector, pcDNA3.2/V5-DEST. These Wnts express the V5 tag.