


This gel depicts the NEB 1kb ladder and a digestion of the plasmid 44167: pEIAV-SIN6.1 CBLacZ-LucWm1A with the following enzymes: HindIII-HF+EcoRI-HF

Plasmid 44167
pEIAV-SIN6.1 CBLacZ-LucWm1A



DNA digested with HindIII-HF+EcoRI-HF double digest for **17 hours**
Gel: 0.8% TAE stained with 100 ml 1X TAE and 5uL of ethidium bromide
Ran for 4 hours at **47 volts**