


Plasmid Name: **pCGN-S2P** ^{WT SS}
 Made by: Jingshi Shen
 Date: 2/1/04

(HA-tagged)



confirmed by sequencing

Method:

vector pCGN XbaI & BamHI

insert part < XbaI S2P 1-163
 BamHI S2P 163-1639 S2P
 template: S2P cDNA clone from library

lig.

For HA ^{SS}(RR) ~~mutant~~ (HA) ^HMA-SSV...
 ↓
 MA-RRV...

Shen et al. 2002
 Dev. Cell 3, 99-111.