

Name	type	location
U6 promoter	promoter	36..192
BsaI	misc_feature	194..199,207..212
sgRNA	misc_RNA	214..293
U6 3'UTR	3'UTR	306..1192
T7 promoter	promoter	rev:1228..1246
M13 fwd	primer_bind	rev:1253..1269
ori	rep_origin	rev:1503..2091
KanR	CDS	rev:2262..3059
bla	promoter	rev:3060..3164
lac promoter	promoter	3497..3527
lac operator	protein_bind	3535..3551
M13 rev	primer_bind	3559..3575
T3 promoter	promoter	3596..3614
SK primer	primer_bind	3651..3654,1..13