


Ple219 Promoter PCR (pEMS1182)**Pleiades Promoter Project**

MiniPromoter: Ple219
pEMS#: 1182
Expected product size (bp): 535

Reaction components	Vol/Rxn (µl)
H ₂ O	15.15
10X PCR buffer*	2.5
50 mM MgCl ₂ *	0.75
2.5 mM dNTPs**	2
10 µM oEMS2646	1.25
10 µM oEMS2622	1.25
Taq Pol. (5 U/µl)*	0.1
DNA***	2
Total Volume of Rxn:	25

L 1 2 3 4 L



L - 100bp Ladder
 1 - Control Plasmid DNA (pEMS1182)
 2 - WT mouse DNA
 3 - No Template
 4 - Transgenic mouse DNA
 Expected band size = 535bp

* Taq Polymerase set from Invitrogen (Cat no.18038-042)

** dNTPs from Invitrogen (Cat no.10297-018)

***Approximately 100 ng DNA used for samples, approximately 5 ng used for plasmid control

Samples run on a 2% agarose gel (containing SYBRsafe (Invitrogen Cat no. S33102))

Cycling conditions:	Step	Temp	Time	Note
	1	94°C	3 min	
	2	94°C	1 min	
	3	61°C	1 min	
	4	72°C	45 sec	repeat steps 2-4 34 times
	5	72°C	5 min	
	6	4°C	hold	

Primers:

Name	Sequence	T _m (°C)	Notes
2646	5'-CTCCTCTTACACACCCAACCTCCCT-3'	56.6	Sense primer for Ple219 in region "3".
2622	5'-GACAGTAGCCACACACCCAGACA-3'	56.3	Antisense primer for Ple219 in region "Prom".