


This gel depicts the NEB 1kb ladder and a digestion of the plasmid N-Myc-53BP1 WT pLPC-Puro with the following enzymes: BamHI-HF and BamHI-HF + XhoI double digest.

**Plasmid 19836**  
N-Myc-53BP1 WT pLPC-



DNA digested with BamHI-HF and BamHI-HF + XhoI double digest overnight  
Gel: 0.8% TAE stained with 100 ml 1X TAE and 5uL of ethidium bromide  
Ran for 2 hours at **60** volts