

**Plasmid Name:** pcDNA-HA-hTRAF3 1-423

**Constructed by:** Gongxian Liao

**Date of Construction:** 10/7/02

**Seletion:** Bacteria: Amp; Mammalian cell: Neo

**Storage:** TRAF3 box


**Characterization:** IB

**Clone # characterized:**

**Cloning Information:**

Human TRAF3 1-423 cDNA was amplified by PCR using SCSH1072/1563 primers, cut with EcoRI/XhoI, and clone into the same sites in pcDNA-HA

**Restriction analysis:**



5' ACC GAT ATC ATG GCC TAC CCT TAT GAT GTG CCA GAT  
KpnI EcoRV HA Tag

TAT GCC TCT CCC GAA TTC ..... CTC GAG AAG CTT TGG  
EcoRI TRAF3 XhoI HindIII

ACT TCT TAA TCT AGA GGG CCC TTC GAA GGT AAG CCT .....3'  
Stop XbaI ApaI BstBI V5 tag