



Comments for pcDNATM 3.3-TOPO[®] 5407 nucleotides

CMV promoter:	47-726
CMV forward primer binding site:	584-604
TOPO [®] cloning site:	741
TK pA reverse primer binding site:	787-805
TK polyadenylation signal:	780-1051
f1 replication origin:	1087-1515
SV40 early promoter:	1520-1889
Neomycin Resistance gene:	1925-2719
SV40 polyadenylation signal:	2895-3025
pUC origin (c):	3408-4081
Ampicillin (<i>bla</i>) resistance gene (c):	4226-5086
<i>bla</i> promoter (c):	5087-5185

(c) = complementary strand