


Comments for pDisplay:  
5325 nucleotides



CMV promoter: bases 1-596  
T7 promoter: bases 638-657  
Murine Ig kappa-chain V-J2-C signal peptide: bases 737-799  
Hemagglutinin A epitope: bases 800-826  
Multiple Cloning Site: bases 827-873  
*myc* epitope: bases 874-903  
PDGFR transmembrane domain: bases 907-1056  
Bovine growth hormone polyadenylation signal: bases 1069-1288  
Col E1 origin: bases 1378-2051  
Thymidine kinase polyadenylation site: bases 2458-2187  
Neomycin/Kanamycin resistance gene: bases 3421-2366  
SV40 origin and promoter: bases 3797-3456  
Ampicillin resistance gene: bases 4736-3876  
f1 origin: bases 5000-5099



Comments for pDisplay™:  
5325 nucleotides