

## pMXs-puro HA-Atg14 dCC\* (siRNA resistance)

Plasmid name: pMXs-puro HA-Atg14 dCC \* (siRNA resistance)

Inserted gene: HA-hAtg14 dCC (lacks the coiled-coil domain (71-184 amino acid residues).

Insert : 1.1kbp

Insert site:EcoRI

Vector: pMXs-puro HA

Antibiotics:Amp

Constructed by: E. Itakura

Reference: MBC Itakura 2008 19, 5360

Human Atg14 siRNA resistant  
siAtg14-1 seq. (582-591)

WT	a gca aat ctt cga cga tcc cat ata ta
siRNA Resistant	a gca aat ctt cgt cgt tca cac ata ta
	A N L R R S H I

