


Ple46 Promoter PCR (pEMS1151)

MiniPromoter: Ple46
pEMS#: 1151
Expected product size (bp): 447

Reaction components	Vol/Rxn (µl)
H ₂ O	15.15
10X PCR buffer*	2.5
50 mM MgCl ₂ *	0.75
2.5 mM dNTPs**	2
10 µM oEMS3043	1.25
10 µM oEMS3037	1.25
Taq Pol. (5 U/µl)*	0.1
DNA***	2
Total Volume of Rxn:	25

L 1 2 3 4 L



L - 100 bp Ladder
1 - Control Plasmid DNA (pEMS1151)
2- WT mouse DNA
3 - No Template
4- Transgenic Mouse DNA
Expected band size = 447 bp

* Taq Polymerase set from Invitrogen (Cat no.18038-042)

** dNTPs from Invitrogen (Cat no.10297-018)

***Approximately 100 ng DNA used for samples, approximately 5 ng used for plasmid control

Samples run on a 2% agarose gel (containing SYBRsafe (Invitrogen Cat no. S33102))

Cycling conditions:	Step	Temp	Time	Note
	1	94°C	3 min	
	2	94°C	1 min	
	3	61°C	1 min	
	4	72°C	45 sec	repeat steps 2-4 34 times
	5	72°C	5 min	
	6	4°C	hold	

Primers:

Name	Sequence	T _m (°C)	Notes
3037	5'- TCAAAGTGCACTCCCAGGCA -3'	56.4	Anti-Sense primer for Ple46 in region "Prom".
3043	5'- CCATCATGGAAAGGGTGTAGTGG -3'	57.0	Sense primer for Ple46 in region "5".