

E-A1 (1 extra tube)


S-A1

Sokol Lab Plasmid Database

| | | | |
|-----------------|--|-------------|-----|
| Plasmid Name | XActB, activinB, bB-64T, bB | Vector Size | |
| Plasmid Size | 4.4 kb | Insert Size | 1.2 |
| Vector Name | pSP64T | mRNA Size | 1.2 |
| Sense/Antisense | Linearize with EcoR1 and transcribe with SP6 for sense XActivin bB mRNA. | | |

V - vectors W- Wnt related
M - markers/probes D- Dsh and below
P - prokaryotic expression constructs
E - eukaryotic expression constructs
O - other constructs

ID E28



Constructed By Sergei Sokol Construction Date 1990
Comments 1.2 kb Nco1 insert patially blunted and cloned into partially blunted Bgl2 site of pSP64T.

Reference Thomsen et al., Cell 1991; Sokol et al., Cell, 1991

Insert Sequence