

Yamanaka Lab's Plasmid Stock # 16-72

Name: pMXs-hLIN28

Size: 5430 bp

Constructed by: K. Takahashi

Date: Nov, 2007

Insert

Description: pENTR-hLIN28

Size: 630 bp

Enzyme used: —

NotI

Human LIN28 ORF


NotI

Vector

Description: pMXs-gw

Size: 4800 bp

Enzyme used: —



Note

NotI digestion

4700, 730 bp

Confirmed by sequencing

(pMXs-51811)

Host: Stb13⁺

Growth Strain

Cloning gateway