


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid **pC0049-EF1a-dPSPCas13b-NES-HIV, H133A/H1058A** with the following enzyme(s): NotI-HF, NheI-HF, KpnI-HF, and SpeI-HF

Plasmid 103865
pC0049-EF1a-dPSPCas13b-NES-HIV,
H133A/H1058A



DNA digested with **NotI-HF/NheI-HF** and **KpnI-HF/SpeI-HF** double digests for **2 hours**.
Gel: **0.8%** agarose stained with 5 uL of SYBR Safe.
Ran for **1.5** hours at **100** volts.