


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pAAV-MBP-2xNLS-tdTomato with the following enzymes: EcoRI-HF and BamHI-HF.

Plasmid 104054:
pAAV-MBP-2xNLS-tdTomato



DNA digested with EcoRI-HF single digest and EcoRI-HF/ BamHI-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 6 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.