


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid pAAV.CMV.PI.EGFP.WPRE.bGH with the following enzymes: HindIII-HF, NotI-HF, and BamHI-HF

**Plasmid 105530**  
pAAV.CMV.PI.EGFP.WPRE.bGH



DNA digested with BamHI-HF single digest and HindIII-HF/NotI-HF double digests for **2 hours**.  
Gel: 0.8% agarose stained with 5 uL of ethidium bromide.  
Ran for 1.5 hours at 100 volts.