


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pAAV.CAG.SF-iGluSnFR.A184V with the following enzymes: SacII and Sall-HF.

**Plasmid 106199:**  
pAAV.CAG.SF-iGluSnFR.A184V



DNA digested with SacII and Sall-HF single digests and SacII/ Sall-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 6 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.