


This gel depicts the NEB 1kb extended ladder and a digestion of the plasmid SunTag-FWAgRNA-4-14aa-TET1cd with the following enzymes: PmeI, SnaBI, HpaI, KpnI-HF, SpeI-HF, and SacI-HF.

Plasmid 106436:
SunTag-FWAgRNA-4-14aa-TET1cd



DNA digested with SpeI-HF and SacI-HF single digests and PmeI/ SnaBI and HpaI/ KpnI-HF double digests for **2 hours**.

Gel: 0.8% agarose stained with 6 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.