

[GENE ID: 26381 Esrrg](#) | estrogen-related receptor gamma [Mus musculus]  
(Over 10 PubMed links)

Score = 484 bits (1246), Expect = 3e-135, Method: Compositional matrix adjust.

Identities = 235/241 (97%), Positives = 237/241 (98%), Gaps = 0/241 (0%)

mutB2	61	LADR <b>A</b> LVVIIGWAKHIPGFSTLSLADQM <del>SL</del> LQSAWMEILILGVVY <b>A</b> SLSFEDEL <del>VY</del> ADDY	120
match		LADR LVVIIGWAKHIPGFSTLSLADQM <del>SL</del> LQSAWMEILILGVVY SLSFEDEL <del>VY</del> ADDY	
ERR3	271	LADR <b>E</b> LVVIIGWAKHIPGFSTLSLADQM <del>SL</del> LQSAWMEILILGVVY <b>R</b> SLSFEDEL <del>VY</del> ADDY	330
mutB2	121	IMDEDQSKLAGLLDLNNAIQLVKKYKSMKLEKEEFVTLKAIALANSDSMHIEDVEAVQK	180
match		IMDEDQSKLAGLLDLNNAIQLVKKYKSMKLEKEEFVTLKAIALANSDSMHIEDVEAVQK	
ERR3	331	IMDEDQSKLAGLLDLNNAIQLVKKYKSMKLEKEEFVTLKAIALANSDSMHIEDVEAVQK	390
mutB2	181	LQDVLHEALQDYEAGQHMEDPRRAGKMLMTLPLL <del>RQT</del> STKAVQHFYNIKLEGKVP <del>MH</del> KLF	240
match		LQDVLHEALQDYEAGQHMEDPRRAGKMLMTLPLL <del>RQT</del> STKAVQHFYNIKLEGKVP <del>MH</del> KLF	
ERR3	391	LQDVLHEALQDYEAGQHMEDPRRAGKMLMTLPLL <del>RQT</del> STKAVQHFYNIKLEGKVP <del>MH</del> KLF	450