


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid pEGFP mCherry-ApoB-EGFP with the following enzymes: HindIII-HF and Agel-HF

Plasmid 112859:
pEGFP mCherry-ApoB-EGFP



DNA digested with: HindIII-HF and Agel-HF single digests and : HindIII-HF/Agel-HF double digests for **2 hours**.

Gel: 0.8% agarose stained with 5 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.