


Sequence: peblindv2.gb (Circular / 4911 bp)
Enzymes: Unique 6+ Cutters (43 of 653 total)
Features: 12 total

Unique Cutters **Bold**



DNA Type: Synthetic DNA

Laboratory Host:

Bacterial Transformation Strain: Unspecified
Dam⁺ Dcm⁺ EcoKI⁺

Description:

Created: Jun 21, 2018

Last Modified: Jun 21, 2018

Accession Number:

Code Number:

Sequence Author:

Comments:

References:

Embedded Files: