


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid pAAV-hSyn-NES-NIR-GECO1 with the following enzymes: Scal-HF and EcoRI-HF

**Plasmid 113683**  
pAAV-hSyn-NES-NIR-GECO1



DNA digested with Scal-HF single digests and Scal-HF/ EcoRI-HF double digests for **2 hours**.  
Gel: 0.8% agarose stained with 5 uL of ethidium bromide.  
Ran for 1.5 hours at 100 volts.