


This gel depicts the NEB 1kb ladder and a digestion of the plasmid AAV-CamKIIa-ChR2-EGFP with the following enzymes: EcoRI-HF and HindIII-HF.

**Plasmid 114216:**  
AAV-CamKIIa-ChR2-EGFP



DNA digested with EcoRI-HF single digest and EcoRI-HF/HindIII-HF double digest for **2 hours.**

Gel: 0.8% agarose stained with 6 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.