

**MASSAGUE LAB PLASMIDS**

CATEGORY: *UCLA Lentiviral*

DATE: *8.12.2018*

PLASMID NAME: ~~TetON-shYAP~~ *TetON-shYAP*

CONSTRUCTED BY: *E.E.E.*

DESCRIPTION: *inducible YAP knock down*

VECTOR: SIZE: kb

INSERT: *YAP shRNA* SIZE: *7* kb *total*

EXPRESSION: *Dox induced* MARKERS: Mammalian: *Puro* Bacterial: Amp

CLONING SITES: *N/A* STORAGE: DNA Prep

CLONING NOTES / SOURCE / REFERENCES:

