


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pCAG PSAM4-5HT3 HC IRES EGFP with the following enzymes: NheI-HF and SacI-HF.

**Plasmid 119740:**  
pCAG PSAM4-5HT3 HC IRES EGFP



DNA digested with NheI-HF single digest and NheI-HF/ SacI-HF double digest for **2 hours**.

Gel: 0.8% agarose stained with 6 uL of ethidium bromide.

Ran for 1.5 hours at 100 volts.