


This gel depicts the NEB 1kb ladder and a digestion of the plasmid 12540 with the following enzymes: BamHI-HF single digest, NotI-HF single digest, and BamHI-HF and NotI-HF double digest

Plasmid 12540

pLPC-N MYC



DNA digested with BamHI-HF single digest, NotI-HF single digest, and BamHI-HF and NotI-HF double digest

2 hours

Gel: 0.8% TAE with 3 uL of ethidium bromide in the gel.

Ran for 4.5 hours at **50** volts