


This gel depicts the NEB 1kb Ladder and a digestion of the plasmid **pEGFP-N1/PD-L1_ECD** with the following enzyme(s): **EcoRI-EF** and **Xho I**

Plasmid 121479
pEGFP-N1/PD-L1_ECD



DNA digested with **EcoRI-HF** single digest and **EcoRI-HF/ XhoI** double digest for **1 hour**.
Gel: **0.8%** agarose stained with **6 uL** of SYBR Safe.
Ran for **1.15** hours at **100** volts.