


This gel depicts the NEB 1kb ladder and a digestion of the plasmid pLX208-gLuc sp-NanoLuc-15aa linker-ICAM-Nrxn3b-HA with the following enzymes: NotI-HF and SacII.

Plasmid 125233:
pLX208-gLuc sp-NanoLuc-15aa linker-
ICAM-Nrxn3b-HA



DNA digested with NotI-HF single digest and NotI-HF/ SacII double digest for **2 hours**.
Gel: 0.8% agarose stained with 6 uL of ethidium bromide.
Ran for 1.5 hours at 100 volts.