


Plasmid name : pMSCVhyg H3.3-WT



Fragment Insert (indicate source, restriction ends, whether sequenced, if PCR derived) : _____

PCR : H3.3-XhoI(5') + H3.3-HpaI(3'); pQCXIPpuro-H3.3-WT-Flag-HA as a template

Antibiotic resistance : Amp _____

Note : Double digestion; XhoI/HpaI (7kb+0.5kb)
 Confirmed by sequencing and Western blot

References : _____

Constructed by Young Hoon Jang _____ Date : 8-6-2015